Duplica^{RealTime} GC/CT/TV/MG Kit
REV: EBR161100_IFU_REV.01D_ENITA

REF: EBR161100- 100 tests

Instructions For Use

INTENDED USE

Duplica^{RealTime} GC/CT/TV/MG Kit is a multiplex Real Time PCR test for the qualitative detection of *Neisseria gonorrhoeae* (GC), *Chlamydia trachomatis* (CT), *Trichomonas vaginalis* (TV) and *Mycoplasma genitalium* (MG) in the urogenital swabs, urine, prostatic liquid and other biological materials.

INTRODUCTION

STDs (Sexually Transmitted Diseases) refer to a variety of bacterial, viral and parasitic infections that are acquired through sexual activity. Some STDs, such as syphilis and gonorrhea, have been known for centuries while others, such as HIV, have been identified only in the past few decades. STDs are caused by more than 25 infectious organisms. As more organisms are identified, the number of STDs continues to expand. Common STDs include: chlamydia, gonorrhea, herpes, HIV, HPV, syphilis, mycoplasma, gardnerella and trichomoniasis.

The development of tests based on nucleic acid amplification technology has been the most important advance in the field of STD diagnosis. Because nucleic acid amplification is exquisitely sensitive and highly specific, it offers the opportunity to use noninvasive sampling techniques to screen for infections in asymptomatic individuals who would not ordinarily seek clinical care.

PRINCIPLE OF THE TEST

Bacterial DNA is extracted from the specimens, amplified and detected by Real Time PCR using fluorescent reporter dye probes specific for *N.gonorrhoeae*, *C.trachomatis*, *T.vaginalis* and *M.genitalium*.

Duplica^{RealTime} GC/CT/TV/MG Kit is a qualitative test that takes advantage of Hot Start PCR, in order to reduce frequency of non specifically primed reaction, and contains the Extraction Amplification Control (EAC). EAC must be used in the extraction procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

REAGENTS PROVIDED

Reagent	Box	Storage (Range, °C)	Volume (ml)	Quantity (tubes)
Amplification Mix (AM)	A	-30 ÷ -16	0,6	1
Taq polymerase (TAQ)	A	-30 ÷ -16	0,06	1
Oligo Mix (OM)*	B	2 ÷ 8	1,1	1
Positive Control (C+)	B	2 ÷ 8	0,2	1
Reaction Blank (BM)	B	2 ÷ 8	0,5	1
Extraction Negative Control (ENC)**	B	2 ÷ 8	1,2	1
Extraction Amplification Control (EAC)***	B	2 ÷ 8	1,0	1

* **store the tube far from light**

** *must be used in the extraction procedure as Negative Control of Extraction and treated as a sample.*

*** *add 10µl of Extraction Amplification Control (EAC) during the DNA extraction procedure directly to the sample/lysis buffer mixture (see herein at the paragraph "DNA purification" and refer to specific kit's protocol).*

STORAGE AND HANDLING

Duplica^{RealTime} GC/CT/TV/MG Kit can be shipped at **2 ÷ 8°C**.

Duplica^{RealTime} GC/CT/TV/MG Kit components are to be stored as follow:

BOX A (AM, TAQ) at **-30 ÷ -16 °C**

BOX B (OM, ENC, BM, EAC, C+) at **2 ÷ 8°C**

The shelf life of reagents before and after the first use is the same, unless otherwise stated.

Duplica^{RealTime} GC/CT/TV/MG Kit is stable up to the expiry date indicated on the kit label. Exposure to light, heat or humidity may affect the shelf life of some components of the kit and should be avoided. Avoid repeated thawing and freezing, as this may reduce the sensitivity of the kit.

MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to specific handbook's section)
- Transport medium
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 µl, 20-200 µl, 100-1000 µl)
- Disposable DNase/RNase-free tips with aerosol barriers
- Tube racks
- Vortex mixer/desktop centrifuge
- PCR box
- Disposable optical polypropylene tubes for Real Time PCR
- Refrigerator
- Deep-freezer
- Thermal cycler for Real Time PCR Rotor-Gene[®] 6000/Q (Corbett Research, Qiagen)

All the equipments should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure optimal performance.

QUALITY CONTROL

In accordance with EuroClone Diagnostica's ISO-9001 and ISO-13485-Certified Quality Management System, each lot is tested against predetermined specifications to ensure consistent product quality.

PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood
- If required, Euroclone Diagnostica offers the necessary technical support for the correct use of the kit
- Carefully read this Instruction for Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipments. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite
- Use powder-free gloves. Do not leave fingerprints on optical part of the tubes
- Materials containing or potentially-containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood
- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet of the product (MSDS, available on request)
- The kit reagents, individual protective equipments, used materials, biological samples and test residuals must be disposed of according to local regulations
- Patient Drug treatment may interfere with the final result of the molecular biology analysis

OPERATING PROCEDURE

a) Sample collection, storage and transport

Sampling of biological materials for PCR analysis, transportation, and storage are described in details in the handbook of the manufacturer. It is recommended that this handbook is read before beginning of the work.

Duplica^{RealTime} GC/CT/TV/MG Kit can analyze DNA extracted from:

- *cervical, conjunctival and urethral swabs*: insert the swab into the DNase/RNase-free 1,5 ml tube and add 0,2 ml of Transport medium. Vigorously agitate swabs in medium for 15-20 sec;
- *urine sediment*: collect 10-20 ml of first-catch urine in a sterile container. Centrifuge for 30 min at 3000 x g, carefully discard the supernatant and leave about 200 µl of solution. Resuspend the sediment. Use the suspension for the DNA extraction;

- *prostatic liquid* stored in "Eppendorf" tube;
- *seminal liquid*: transfer about 30 µl of seminal liquid to a polypropylene tube (1,5 ml) and add 70 µl of sterile saline solution.

It is recommended to process samples immediately after collection. Store samples at 2 ÷ 8°C for no longer than 24 hours, or freeze at -30 ÷ -16 C. Transportation of clinical specimens must comply with country, federal, state and local regulations for the transport of etiologic agents.

b) DNA purification

EuroClone Diagnostica recommends to use Duplicα EasyBind DNA1 kit (ref. EDR005050) for manual purification. Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

Add **10 µl** of **Extraction Amplification Control** (EAC) during the DNA isolation procedure directly to the sample/lysis buffer mixture.

c) Thermal Cycler setup

Refer to the specific handbook of the equipment used to set the Thermal Profile.

We recommend to switch on the instrument, and to set the thermal profile before preparing the reaction mix.

Thermal Profile

Step	Temperature, °C	Time	Repeats
Hold	95	15 min	1
Cycling 1	95	5 s	5
	60	20 s	
	72	15 s	
Cycling 2	95	5 s	40
	60	20 s <i>Fluorescence acquisition</i>	
	72	15 s	

Fluorescence is detected at the second step of Cycling 2 stage (60°C) in Green, Yellow, Orange, Crimson and Red fluorescence channels

d) Preparation of PCR mix

The total reaction volume is **25 µl**.

Thaw Amplification Mix (AM) before mixing.

For each experiment prepare a Reaction Mix considering 2 controls (**ENC** and **C+**), 1 Reaction Blank (**BM**) and **n+1** samples.

The reagents have to be mixed as indicated in the table below:

REAGENT	VOLUME (µl)
Oligo Mix (OM)	10
Amplification Mix (AM)	5
Taq Polymerase (TAQ)	0,5

The PCR Mix has to be freshly prepared every time

Vortex and centrifuge for 2-3 sec.

After its preparation, aliquot **15 µl** of **Reaction Mix** in the tubes, then add **10 µl** of **extracted DNA** (or **ENC, BM, C+**) sample to appropriate tube. Mix by pipetting. Place the tubes in the instrument and start the program of amplification.

e) Instrument settings

Rotor-Gene® 6000/Q

Make the adjustment of the fluorescence channel sensitivity:

1. *Channel Setup* → *Gain Optimisation* → *Auto Gain Optimisation Setup* → *Optimise Acquiring* and select *Perform Optimisation Before 1st Acquisition*
2. For *Green* channel indicate *Min Reading 5, Max Reading 10*
3. For *Yellow, Orange, Red, Crimson* channels *Min Reading 4, Max Reading 8*
4. In the column *Tube position* program position of the tubes in the carousel of the Rotor-Gene
The first position must contains reaction tube with reagents
5. Close the window *Auto Gain Calibration Setup*

f) ANALYSIS and INTERPRETATION of RESULTS

- **Neisseria gonorrhoeae** DNA is detected on the **Green** channel,
- **Chlamydia trachomatis** DNA is detected on the **Yellow** channel,
- **Mycoplasma genitalium** DNA is detected on the **Orange** channel,
- **Trichomonas vaginalis** DNA is detected on the **Crimson** channel,
- **Extraction Amplification Control (EAC)** is detected on the **Red** channel.

The results are interpreted with the software of **Rotor-Gene®** through the presence of crossing of fluorescence curve with the threshold line:

1. Press *Analysis* then select button *Quantitation*.
2. Perform the operation for the channel **Green** (*Cycling A.Green*), then refer to the **Settings Table** for the channels **Yellow** (*Cycling A.Yellow*), **Orange** (*Cycling A.Orange*), **Red** (*Cycling A.Red*) and **Crimson** (*Cycling A.Crimson*)

2.1. Data analysis of *Neisseria gonorrhoeae* DNA

- Select **Green** channel window.
- Select the **Dynamic tube** in the main window menu.
- In **CT Calculation** menu set **Threshold = 0,1**.
- Select **Outlier Removal** button and type **0** in the text field.
- The C_T (Threshold cycle) values for each sample in channel will be shown in the results grid.

Settings Table

	<i>C.trachomatis</i> DNA	<i>M.genitalium</i> DNA	<i>T.vaginalis</i> DNA	Extraction/Amplification Control (EAC)
Channel Window	<i>Yellow</i>	<i>Orange</i>	<i>Crimson</i>	<i>Red</i>
Dynamic Tube	<i>On</i>	<i>On</i>	<i>On</i>	<i>On</i>
Threshold in CT calculation	<i>0,1</i>	<i>0,1</i>	<i>0,1</i>	<i>0,07</i>
Outlier Removal	<i>5</i>	<i>5</i>	<i>10</i>	<i>5</i>

The results of the analysis are reliable only if the results obtained for both Positive and Negative Controls are correct. Refer to the here below **Run Validation** table.

Run Validation

Control	Stage for control	C_T Green	C_T Yellow	C_T Orange	C_T Crimson	C_T Red	Interpretation
ENC	DNA isolation	Neg	Neg	Neg	Neg	<33	Valid result
BM	Amplification	Neg	Neg	Neg	Neg	Neg	Valid result
C+	Amplification	<35	<35	<35	<35	<33	Valid result

If the above conditions have been met the run is valid and it is possible to analyze the data, refer to following instructions:

1. The sample is considered to be **positive** for ***Neisseria gonorrhoeae*** if its C_T value is defined (the fluorescence curve crosses the threshold line) in the **Green** channel.
2. The sample is considered to be **positive** for ***Chlamydia trachomatis*** if its C_T value is defined (the fluorescence curve crosses the threshold line) in the **Yellow** channel.
3. The sample is considered to be **positive** for ***Mycoplasma genitalium*** if its C_T value is defined (the fluorescence curve crosses the threshold line) in the **Orange** channel.
4. The sample is considered to be **positive** for ***Trichomonas vaginalis*** if its C_T value is defined (the fluorescence curve crosses the threshold line) in the **Crimson** channel.
5. The sample is considered to be **negative** for ***Chlamydia trachomatis*, *Neisseria gonorrhoeae*, *Mycoplasma genitalium* and *Trichomonas vaginalis*** if its C_T value is not defined (the fluorescence curve does not cross the threshold line) in Green, Yellow, Orange and Crimson channels and the C_T value does not exceed the boundary value in the Red channel ($C_t < 33$).

QUALITY CONTROL PROCEDURE

A defined quantity of Extraction Amplification Control (EAC) is introduced at the beginning of sample preparation procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

Extraction Negative Control (ENC), Reaction Blank (BM), Positive Control (C+) are required for every run to verify that the specimen preparation, the amplification and the detection steps are performed correctly.

If the controls are out of their expected range (see **Run Validation** table), all specimens and controls from that run must be processed again beginning from the sample preparation step.

PERFORMANCE

Sensitivity

The analytical sensitivity for *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, and *Trichomonas vaginalis* DNA is not less than 5×10^2 Genome Equivalents per 1 ml of sample (GE/ml).

The analytical sensitivity of each microorganism does not change even if three other microorganisms are present at high concentrations.

Specificity

The analytical specificity of **Duplica^{RealTime} GC/CT/TV/MG Kit** is ensured by selection of specific primers and probes as well as stringent reaction conditions. The primers and probes were checked for possible homologies to all sequences published in gene banks by sequence comparison analysis. The clinical specificity of **Duplica^{RealTime} GC/CT/TV/MG Kit** was confirmed in laboratory clinical trials.

TROUBLESHOOTING

Problem 1: Weak or no signal of the Extraction Amplification Control (EAC) in the Red channel.

1. The PCR was inhibited:
 - Make sure that you use a recommended DNA extraction method and follow to the manufacturer's instructions
 - Centrifuge for 2 min at maximum speed the extracted DNA and take carefully the supernatant for the PCR; sorbent present in the pellet inhibits the PCR reaction.
2. The reagents storage conditions didn't comply with the instructions:
 - Check the storage conditions.
3. Improper DNA extraction:
 - Repeat analysis starting from the DNA extraction stage.
4. The PCR conditions didn't comply with the instructions:
 - Check the PCR conditions and select for the the Extraction Amplification Control detection the fluorescence channel reported in the protocol.
5. The Extraction Amplification Control was not added to the sample during the pipetting of reagents:
 - Make attention during the DNA extraction procedure.

Problem 2: Weak or no signal of the Positive Control (C+).

1. The PCR conditions didn't comply with the instructions:
 - Check the amplification protocol and select the fluorescence channel reported in the manual.

Problem 3: Any signal with Extraction Negative Control (ENC) except for Red channel.

1. Contamination during DNA extraction procedure. All samples results are INVALID:
 - Decontaminate all surfaces and instruments with sodium hypochlorite (0,5%) and ethanol (70%) or special DNA decontamination reagents
 - Use only filter tips during the extraction procedure. Change tips between tubes
 - Repeat the DNA extraction with the new set of reagents.

Problem 4: Any signal with Reaction Blank (BM).

1. Contamination during PCR preparation procedure. All samples results are INVALID:
 - Decontaminate all surfaces and instruments with sodium hypochlorite (0,5%) and ethanol (70%) or special DNA decontamination reagents
 - Pipette the Positive control at last
 - Repeat the PCR preparation with the new set of reagents.

Problem 5: Fluorescence intensity varies.

1. The PCR Master Mix is not well prepared:
 - Carefully repeat the PCR preparation procedure
2. Air bubbles trapped in the PCR tubes:
 - Check the presence of air bubbles before starting a new run.

Problem 6: Absence of any fluorescent signal.

1. Verify the performance of the thermalcycler:
 - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions didn't comply with the instructions:
 - Check the storage conditions

- Check the expiry date of the kit.

Problem 7: The thermalcycler gives an error message.

1. Refer to the Real Time PCR instrument user manual or contact the local technical support.

REFERENCES (see pag. 12/12)

"THE PURCHASE OF THIS PRODUCT GRANTS THE PURCHASER RIGHTS UNDER CERTAIN ROCHE PATENTS TO USE IT SOLELY FOR PROVIDING HUMAN IN VITRO DIAGNOSTIC SERVICES. NO GENERAL PATENT OR OTHER LICENSE OF ANY KIND OTHER THAN THIS SPECIFIC RIGHT OF USE FROM PURCHASE IS GRANTED HEREBY"

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EuroClone Diagnostica Srl