DEUTSCH

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REFERENCES

1. Lake IR, *et al.*, Case-control study of environmental and social factors influencing cryptosporidiosis. *Eur J Epidemiol.* 2007; 22:805-811
2. Khalakdina A, *et al.*, Is drinking water a risk factor for endemic cryptosporidiosis? A case-control study in the immunocompetent general population of the San Francisco Bay Area. *BMC Public Health.* 2003; 3:1471-2458
3. Cacciò SM, *et al.*, Unravelling *Cryptosporidium* and *Giardia* epidemiology. *Trends Parasitol.* 2005; 21:430-437
4. Certad G, *et al.*, *Cryptosporidium parvum*, a potential cause of colic adenocarcinoma. *Infect. Agent. Cancer* 2007; 2:22
5. Nygard K, *et al.*, A large community outbreak of waterborne giardiasis-delayed detection in a non-endemic urban area. *BMC Public Health.* 2006; 6: 141
6. Monis PT, Thompson RCA. *Cryptosporidium* and *Giardia*-zoonoses: fact or fiction? *Infect Genet. Evol.* 2003; 3:233-244
7. Ekdahl K, Andersson Y. Imported giardiasis: impact of international travel, immigration, and adoption. *Am. J. Trop. Med.* 2005; 72: 825-830

8. Ungar BLP, *et al.*, Use of a monoclonal antibody in an enzyme immunoassay for the detection of *Entamoeba histolytica* in fecal specimens. Am. J. Trop. Med. Hyg. 1985; 34:465-472
9. Strachan WD, *et al.*, Immunological differentiation of pathogenic and non-pathogenic isolates of *Entamoeba histolytica*. The Lancet 1988; 12:561-563
10. Yau YCW, *et al.*, Development of monoclonal antibodies which specifically recognise *Entamoeba histolytica* in preserved stool samples. 39: 716-719.
11. Braga LL. Seropositivity for intestinal colonization with *Entamoeba histolytica* and *Entamoeba dispar* in individuals in Northeastern Brazil. 36: 3044-3045.
12. Petri Jr WA, *et al.*, Pathogenic and nonpathogenic strains of *Entamoeba histolytica* can be differentiated by monoclonal antibodies to the galactose-specific adherence lectin. 1990; 58:1802-1806
13. Paniagua GL, *et al.*, Two or more enteropathogens are associated with diarrhoea in Mexican children. 2007; 28: 17
14. Haque R, *et al.*, *Entamoeba histolytica* infection in children and protection from subsequent amebiasis. 2006; 74: 904-909

REV. 05/15

 Manufactured By	<p>Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124</p>
 Authorized Representative	<p>Meridian Bioscience Europe S. r. L Via dell' Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu</p>

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleud
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info.bnl@meridianbioscience.eu

Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 411 62 11 66
 Fax: +31 (0) 411 62 48 41
 Email: info.bnl@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		DET REAG	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.