

M. pneumoniae + *C. pneumoniae* Duplex Real-Time PCR



25 / 100 / 150 tests



Instructions for use



Labsystems Diagnostics Oy
Tiilitie 3, FIN-01720 Vantaa, Finland
Tel. +358-20-155 7530, Fax +358-20-155 7521
E-mail: contact@labsystemsdx.com
www.labsystemsdx.com

25.09.2015

Instructions for use

For *in vitro* diagnostic use only

A quantitative Real-Time PCR assay for the rapid detection of *Mycoplasma pneumoniae* and *Chlamydophila pneumoniae* in human respiratory samples.

Products no.:
8100 100, 25 wells
8101 100, 100 wells
8100 101, 150 wells

CONTENTS

	Page
INTENDED USE	2
INTRODUCTION	2
PRINCIPLE OF THE TEST	2
KIT CONTENTS	2
REAGENT PREPARATION	3
MATERIALS REQUIRED BUT NOT PROVIDED	3
PRECAUTIONS	3
SAMPLE COLLECTION AND HANDLING	4
TEST PROCEDURE	4
LIMITATION OF THE PROCEDURE	5
PERFORMANCE CHARACTERISTICS	5
TROUBLE SHOOTING	6
NOTICE TO THE PURCHASER	6
REFERENCES	6
RELATED PRODUCTS	7
SYMBOLS USED	8

INTENDED USE

Labsystems Diagnostics' *M. Pneumoniae + C. pneumoniae* Real-Time PCR kit is developed for the rapid and simultaneous detection of *Mycoplasma pneumoniae* and *Chlamydophila pneumoniae* in human samples. Kit includes an cellular control detecting human genomic DNA to insure the presence of cells in the samples.

The positive result is an aid for the diagnosis of acute *Mycoplasma pneumoniae* and *Chlamydophila pneumoniae* infection.

Kit is for professional use only.

INTRODUCTION

Acute respiratory tract infections are widespread in adults and children. *Chlamydophila pneumoniae* is primarily a respiratory tract pathogen that causes approximately 10-20% of community acquired pneumonia in adults and children, and 10-20% of acute bronchitis in adults (1, 2, 3). It also causes sinusitis, primary pharyngitis, and may trigger for asthma (4). Most of infections with this micro-organism are in fact subclinical and asymptomatic and only rarely cause on overt disease (2). Chronic infection with *C. pneumoniae* has been suggested as a factor in the development of atherosclerosis (5, 6).

Mycoplasma pneumoniae is a common human pathogen causing asymptomatic, mild or, rarely, severe upper and lower respiratory tract infections. *M. pneumoniae* infection has been found to account for 15-20% of all cases of pneumonia (7). The infection with *M. pneumoniae* is usually endemic, however small epidemic peaks are observed at several years' intervals (7). *M. pneumoniae* is weakly contagious affecting mainly children, young adults and immunosuppressed adults (8). Because *M. pneumoniae* lacks the cell wall, commonly prescribed β -lactam antibiotics are not effective. Therefore the correct diagnosis is essential to adapt the treatment to shorten the duration of the disease.

To date, most investigations have relied on serologic immunoassays or culture tests. Those technical problems are avoided with the Real-Time PCR assay developed by Labsystems Diagnostics, providing easy, fast and sensitive results.

PRINCIPLE OF THE TEST

The principle of the Labsystems Diagnostics' *M. Pneumoniae + C. pneumoniae* Real-Time PCR kit is based on the 5' nuclease technique. The kit is composed of ready to use optimized reaction mixture for the amplification and detection. This kit is intended to use on respiratory samples after DNA extraction (not included in the kit).

The kit is suitable for all Real-Time PCR instruments with minimum of three measurement channels (FAM, VIC and ROX). Kit is suitable for both fast instruments (such as ABI 7500FAST Real-Time PCR system from Life Technologies; PikoReal24 from Thermo Scientific; Rotor-Gene Q from Qiagen) as well as standard instruments (such as ABI Prism 7500, 7700 and 7900 from Life Technologies). Kit does not contain passive reference dye.

The target sequences identified by the kit have been selected from the conserved regions of the pathogen of interest.

Target

Mycoplasma pneumoniae
Chlamydophila pneumoniae

Gene

P1 cytadhesine gene
ompA major outer membrane protein gene

The amplification of *Mycoplasma pneumoniae* is measured with the FAM fluorescence channel whereas the amplification of *Chlamydophila pneumoniae* is measured with the VIC/HEX fluorescence channel. The amplification of the internal control is measured with the ROX channel. The assay can be used in any Real-Time PCR instrument provided with the detection channels for the mentioned fluorescent labels.

KIT CONTENTS

Note:

- Reagents are stored at between +2°C and +8°C, in to separate area.
- The expiration date is printed on each component label and on the package. Do not use reagents after the expiration date.
- Protect from light. The light sensitive amplification mix is packaged in brown vials for protection.

- To avoid contamination, store the positive control before and after first opening in the same area as the extracted samples.
- 1 AMPLIFICATION MIX, 1 tube (700 µl)/4 tubes/6 tubes
Contains dNTPs, MgCl₂, PCR buffer, specific primers and probes, polymerase enzyme
- 2 NEGATIVE CONTROL, 1.0 ml
Molecular biology grade sterile water. The control must be amplified at the same time and using the same protocol as the patient samples. The negative control verifies the absence of contamination during amplification.
- 3 POSITIVE CONTROL, 20 µl / 2x 20 µl
Contains two plasmids which are amplified with the primers and probes targeted for *M. and C. Pneumoniae* in the amplification mix. It verifies the proper performance of the primers and probes. Signal is detected with the same parameters as for the targets. Control can be used to generate a standard curve for copy number determination.

Standard curve	Copy Number
Tube 1 Positive Control (Vial 3)	200 000 per µl
Tube 2	20 000 per µl
Tube 3	2 000 per µl
Tube 4	200 per µl
Tube 5	20 per µl
Tube 6	2 per µl

- 6) Pipette 5 µl of standard template into each well, according to your experimental plate set up. The final volume in each well is 30 µl. Do not store and re-use the diluted control.

MATERIALS REQUIRED BUT NOT PROVIDED

- DNA extraction kit
- Real-Time PCR instrument with at least three detection channels suitable for FAM, HEX/VIC and ROX measurement
- Disposable gloves
- Adjustable pipettes
- Sterile pipette tips with filter
- Vortex mixer
- Centrifuge with a rotor for 2.0 ml tubes
- Sterile nuclease-free microtubes for positive control dilution series
- Sterile nuclease-free PCR tubes, strips or plates recommended by the Real-Time PCR instrument manufacturer

REAGENT PREPARATION

Reagent	Preparation	Stability of opened /diluted reagents (+2°C to +8°C)
1 Amplification mix	Ready for use	6 months *)
2 Negative control	Ready for use	12 months *)
3 Positive control	Ready for use; optional for quantitation: dilute the control (vial 3) with Negative control water (vial 2) for the ten-fold dilution series standard curve	6 months *)

- *) The stability of the opened reagents is the maximum only if they are stored properly at +2°C to +8°C. High environmental temperature and contamination may decrease the stability. Do not re-use diluted positive control.

Preparation of standard curve dilution series

- 1) Pipette 90 µl of nuclease free sterile water (Negative Control, Vial 2) into 5 tubes and label 2-6.
- 2) Pipette 10 µl of Positive Control (Vial 3) into tube 2
- 3) Vortex thoroughly
- 4) Change pipette tip and pipette 10 µl from tube 2 into tube 3
- 5) Vortex thoroughly

Repeat steps 4 and 5 to complete the dilution series.

PRECAUTIONS

For *in vitro* diagnostic use only.

Extreme precautions should be taken to prevent contamination of the reactions. Amplification procedures require highly skilled professionals.

Use separate working places for amplification reactions and sample preparations. To avoid contamination, three separate working areas are recommended. Clean area I to store the kit, area II for sample preparation and positive control and designated area III to perform the amplification. Allocate a set of lab coats and pipettes to each area. Never introduce an amplified product in reagent and/or sample preparation areas.

Samples must be prepared under a biological safety hood.

Discard all materials and samples as if capable of transmitting infection. The preferred method of disposal is autoclaving for a minimum of one hour at 121°C. Liquid wastes not containing acid and neutralized waste may be mixed with sodium hypochlorite in volumes such that the final mixture contains 50-500 mg/l free available chlorine. Allow 30 minutes for decontamination to be completed. Spills should be wiped off thoroughly using either an iodophor disinfectant or sodium hypochlorite solution. Materials used to wipe off spills should be added to biohazardous waste matter for proper disposal.

Tubes from specimens and amplification mix must never be opened at the same time.

Wear disposable gloves while handling samples and kit reagents. Afterwards wash hands carefully. Never pipette by mouth.

Do not use same pipettes to handle samples and amplification mix. All pipettes used must be equipped with sterile filter tips.

Do not reuse a strip or plate even if some wells were not used.

Accurate and precise pipetting, as well as following the exact time and temperature requirements, is essential.

Do not eat, drink or smoke in dedicated work areas.

SAMPLE COLLECTION AND HANDLING

Clinical samples should be transported to the laboratory as soon as possible, aliquoted and stored at -70°C or -20°C until use. **Samples should not be repeatedly frozen and thawed.**

Microbially contaminated samples may give erroneous results.

The following methods can be used for DNA extraction of clinical specimens. For manual extractions, use the QIAamp DNA Mini Kit (Qiagen) or QIAamp MinElute Virus Spin Kit (Qiagen) according to the manufacturer's instructions. Briefly, extract 200 µl of sample material and elute the nucleic acid in 50 – 100 µl. Automated DNA extraction can be performed for example with the KingFisher nucleic acid purification system (Thermo Scientific) or MagNA Pure 96 instrument (Roche).

Samples should be processed within one week after the DNA isolation.

Sample types evaluated with this kit: sample transfer medium, nasal aspiration biopsy, respiratory biopsy, nasal swabs and bronchoalveolar lavage.

TEST PROCEDURE

PRELIMINARY PREPARATIONS

- **Prepare the standard curve dilution series before starting the assay.**
- **The use of standard curve is optional. When using single well control, perform the first dilution step (10 µl positive control + 90 µl nuclease free sterile water) and use 5 µl of that dilution as positive control sample.**

Prestep Dilute the samples and controls if necessary to sterile water or negative control (vial 2).

PCR Reaction set-up

The PCR reaction is performed in a final volume of 30 µl. The composition of the amplification reaction is as follows:

Amplification mix	25 µl
Extracted sample or control	5 µl
Total	30 µl

PCR amplification program

After placing the PCR tubes into the instrument the following protocol is used:

TEMPERATURE	TIME	CYCLES
95°C	1 min	1
95°C	3 sec	45
60°C	35 sec	
20°C	10 sec	1

Note 1: Select "none" in "passive reference" or "reference dye" when programming.

Note 2: Detection in the FAM, HEX/VIC and ROX channels should be activated.

Detailed programming guides for each Real-Time PCR instrument are available upon request.

NOTES:

1. Use of duplicates is preferable. **Mix well.**
2. **Avoid contamination:** When removing aliquots from the reagent vials, use aseptic technique to avoid contamination. Use a new pipette tip for each sample. Do not touch the walls of the wells with pipette tips.
3. Perform at least one negative and one positive control amplification in every PCR run to verify the reaction performance.
4. Do not store and handle Amplification mix in the same area as the samples and positive control to avoid contamination.

RESULTS

Programming of the instruments is carried out according to the user's manual. Fluorescence data are plotted against the number of cycles. The threshold cycle (C_t or C_q) serves as a tool for calculation of the starting template amount in each sample. The threshold is adjusted to a value above the baseline, but must be located in the log-linear range of the PCR. Before determining the C_t value, ensure that the baseline is positioned correctly and adjust if necessary.

The negative control should not give a detectable signal on any channel.

The positive control should give a detectable signal at both FAM and VIC channels below 30 cycles.

Interpretation of the Results

A positive sample displays a C_t value.

If a C_t value can not be calculated, the sample is considered as negative or inhibited or contains a limited number of target cells.

No C_t value

$35 \leq C_t \leq 44$

$C_t < 35$

Negative

Low positive

Positive

When a low positive sample is detected, sample should be retested.

		<i>Chlamydomphila pneumoniae</i>	
		Ct value detected	No Ct value
<i>Mycoplasma pneumoniae</i>	Ct value detected	<i>Chlamydomphila pneumoniae</i> + <i>Mycoplasma pneumoniae</i> +	<i>Chlamydomphila pneumoniae</i> – <i>Mycoplasma pneumoniae</i> +
	No Ct value	<i>Chlamydomphila pneumoniae</i> + <i>Mycoplasma pneumoniae</i> –	<i>Chlamydomphila pneumoniae</i> – <i>Mycoplasma pneumoniae</i> –

The quality of the sample taken and/or inhibition for the amplification reaction is verified using the internal control. If the C_t value is less than 32, the DNA extraction is performed correctly and the number of cells used in the assay is adequate.

LIMITATION OF THE PROCEDURE

Because no single method leads to the definitive diagnosis, the results of the present method should be interpreted in conjunction with the clinical condition, epidemiological situation and other laboratory methods.

It is recommended that the assay is performed by qualified and trained laboratory technician.

PERFORMANCE CHARACTERISTICS

Analytical specificity

The specificity of the primers and probes in the *M. pneumoniae* + *C. pneumoniae* Real-Time PCR kit was determined by sequence analysis (viral, bacterial and human). The kit was experimentally evaluated with real-Time PCR on the following pathogens:

Bordetella pertussis

Bordetella parapertussis

Chlamydomphila pneumoniae (strain CM-1 and strain AR-39)

Chlamydia trachomatis (strains UW-36/Cx, UW-43/Cx, BOUR)

Haemophilus influenzae

Influenza A, B

Legionella pneumophila subs. *Pneumophila*

Mycobacterium tuberculosis

Mycoplasma pneumoniae

Mycoplasma bovis

Mycoplasma genitalium

Mycoplasma fermentas

Mycoplasma arginini

Streptococcus pneumoniae

No cross-reaction between *Chlamydomphila pneumoniae* and *Mycoplasma pneumoniae* and no amplification with any of the other pathogens were observed.

Tests were also performed on human DNA extracts negative for *Chlamydomphila pneumoniae* and *Mycoplasma pneumoniae*. No amplification of sequences of human origin were observed.

External evaluation

Clinical performance evaluation of Labsystems Diagnostics' *M. pneumoniae* + *C. pneumoniae* Duplex Real-Time PCR was performed in European Hospital Laboratory. Altogether 92 respiratory patient samples were analyzed with 100 % concordance with the in-house diagnostic Real-Time PCR assay.

n = 91*		Clin lab in-house		
		MP+	CP+	NEG
Labsystems	MP+	25	0	0
	CP+	0	27	0
	NEG	0	0	40

* one sample was tested positive for both MP and CP

Sample content	<i>M. pneumoniae</i> positive (%)	<i>C. pneumoniae</i> positive (%)
<i>M. pneumoniae</i> (24 samples)	24 (100%)	0
<i>C. pneumoniae</i> (26 samples)	0	26 (100%)
<i>M.+C. pneumoniae</i> (1 sample)	1 (100%)	1 (100%)
Negative (40 samples)	0 (100%)	0 (100%)

Quality control panels

Assay specificity and sensitivity is tested regularly using external quality control panels.

TROUBLESHOOTING

Cause/Error	Solution
Positive control is negative	
1. One of the reaction components was not added	Ensure that all components have been added
2. The positive control was not stored properly	Store all components according to manufacturers instructions
3. Positive control has been degraded	Follow the general manufacturers instructions
4. Wrong PCR profile	Check the programming of the Real-Time PCR instrument

Cause/Error	Solution
Negative control is positive	
1. Carry over contamination	Repeat the entire experiment with fresh reagents. Handle samples, kit components and consumables as prescribed.

Cause/Error	Solution
Internal control is negative	
1. No DNA in the sample	Repeat extraction procedure to ensure the DNA extraction. Choose another DNA isolation method.
2. Inhibitory substances	Dilute DNA sample and repeat the experiment. Repeat the extraction procedure to ensure that no contaminants are present.

Cause/Error	Solution
Very weak fluorescence signal	
1. Incorrect instrument settings	Check channel settings
2. Improper storage of the Amplification mix	Store all components according to manufacturers instruction. Protect the Amplification mix from light exposure.

NOTICE TO THE PURCHASER

LIMITED LICENSE.

The purchase of this product includes a license to use only this amount of product solely for the purchaser's own use solely in the human in vitro diagnostic field (in accordance with applicable FDA and other regulatory requirements) and may not be used for any other commercial use, including without limitation repacking or resale in any form. No general patent or other license of any kind other than this specific right of use from purchase is granted hereby by Labsystems Diagnostics Oy.


QIAGEN and QIAamp are registered trademarks of the QIAGEN Group.
FAM, HEX, VIC, ROX, ABI and ABI Prism are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries.
KingFisher and PikoReal are registered trademarks of Thermo Fisher Scientific.
LightCycler is a registered trademark of Roche Group.
Rotor-Gene is a trademark of Corbett Life Science.

REFERENCES:

- Grayston JT. *Chlamydia pneumoniae*, strain TWAR. *Chest* 1989; 95:664-9.
- Kleemola M, Saikku P, Visakorpi R, Wang SP, Grayston JT. Epidemics of pneumonia caused by TWAR, a new Chlamydia organism, in military trainees in Finland. *J Infect Dis* 1988;157(2):230-6.
- Grayston JT, Aldous MB, Easton A, Wang SP, Kuo CC, Campbell LA, Altman J. Evidence that *Chlamydia pneumoniae* causes pneumonia and bronchitis. *J Infect Dis* 1993; 168(5):1231-5.
- Hahn DL, Dodge RW and Golubjatnikov R. Association of *Chlamydia pneumoniae* (strain TWAR) infection with wheezing, asthmatic bronchitis and adult-onset asthma. *JAMA* 1991;266(2):225-30.
- Saikku P, Leinonen M, Mattila K, Ekman MR, Nieminen MS, Mäkelä PH *et al.* Serological evidence of an association of a novel Chlamydia, TWAR, with chronic coronary heart disease and acute myocardial infarction. *Lancet* 1988;2(8618):983-6.
- Saikku P, Leinonen M, Tenkanen L, Linnanmäki E, Ekman MR, Manninen V *et al.* Chronic *Chlamydia pneumoniae* infection as a risk factor for coronary heart disease in the Helsinki Heart Study. *Ann Intern Med* 1992;116(4):273-8.
- Foy HM. Infections caused by *Mycoplasma pneumoniae* and possible carrier state in different populations of patients. *Clin Infect Dis* 1993;17, Suppl.1:S37-46.
- Taylor-Robinson D. Infections due to species of *Mycoplasma* and *Ureaplasma*: An Update. *Clin Infect Dis* 1996;23:671-84.



CE-mark



CE-mark, code of the Notified Body



Catalog number



Contains sufficient for < n > tests



Use by YYYY-MM



Batch code



Temperature limitation



In vitro diagnostic medical device



Manufacturer



Consult instructions for use



Positive control



Negative control



Potential biohazardous material