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Instructions For Use

Seraline[®] Zöliakie-3 IgG Seraline[®] Zöliakie-3 IgA

Line Immunoassay for detection of IgG or IgA antibodies
to Gliadin, Tissue Transglutaminase (tTG) and Mannan in human serum or plasma

 LIA-001-3 G  LIA-001-3 A  20  *In-vitro*-Diagnostic device 

Seramun Diagnostica GmbH · Spreehagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Distribution: Sekisui Virotech GmbH · Löwenplatz 5 · 65428 Rüsselsheim · Germany · www.sekisuivirotech.com
phone +49 (0) 6142 6909 0 · fax +49 (0) 6142 9666 13 · info@sekisuivirotech.com

Introduction

Coeliac disease (also called non-tropical sprue, coeliac sprue, gluten sensitive enteropathy and gluten intolerance) is a digestive disorder in genetically-predisposed individuals. The only susceptibility loci established are the HLA-DQ2 and DQ8. Coeliac disease is characterised by damage or flattening of all or part of the villi lining the small intestine, causing scar tissue that cannot absorb nutrients. This damage is caused by exposure to gluten (gliadin) and related proteins found in wheat, rye, malt, barley and oats. Damage to the villi reduces the ability of the intestines to absorb nutrients, and it is believed that the resulting nutritional deficiencies likely cause the wide spectrum of symptoms associated with the disorder. Coeliac disease may lead to digestive problems, such as indigestion, heartburn and irritable bowel syndrome, unexplained weight loss or other signs of nutritional deficiency due to malabsorption, and a wide range of other problems in different bodily functions. In young children, the most common symptoms are steatorrhoea, weight loss, abdominal distension, and slow growth/failure to thrive, but irritability, vomiting and tiredness are common. It has been suggested that some cases of autism may be caused by coeliac disease. The gold standard test for coeliac disease is still upper endoscopy with biopsy of the distal duodenum or jejunum. To avoid false negative results, the first endoscopy must be done while the patient is on a normal, gluten-containing diet or very shortly after going on a gluten-free diet. Serology has been proposed as a screening measure, because the presence of serum IgA antibodies reactive against Gliadin and Tissue Transglutaminase (tTG) is indicative of coeliac disease. Like endoscopy, these tests are not accurate in patients who have been on a gluten-free diet for some time; they must be performed while the person is on a normal diet or within a few months after eliminating gluten from the diet.

In children younger than 2 years suffering from coeliac disease antibodies to Tissue Transglutaminase are often not detectable. Therefore combined testing for anti-Gliadin and anti tTG-antibodies is recommendable. In adults, the symptoms of coeliac disease may be mistaken for irritable bowel syndrome (IBS) or an inflammatory bowel disease such as Crohn's disease. The latter can be diagnosed with a sensitivity up to 75 % by detection of anti- *Saccharomyces cerevisiae* antibodies (ASCA).

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Intended use

The *Seraline*[®] Zöliakie-3 IgG / *Seraline*[®] Zöliakie-3 IgA is a Line Immunoassay (LIA) for detection of specific IgG or IgA antibodies against the following antigens: Gliadin, Tissue Transglutaminase (tTG) and Mannan in human serum or plasma samples.

Principle of the test

The antigens are plotted on nitrocellulose membranes and blocked to avoid unspecific reactions. After this procedure the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

Step 1

The test strips are incubated in the diluted (1 : 101) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates.

The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template. Alternatively the strips can be evaluated with the *Seraline*[®]scan software.

Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 101 (15 µl sample and 1500 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

Kit components

1	TESTSTR	test strips incl. evaluation template	20 strips colour coded: yellow
2	WIB CONC 5x	wash and incubation buffer 5-fold	70 ml concentrate, transparent bottle, black cap
3	CONJ HRP IgG CONJ HRP IgA	anti-human IgG-HRP-conjugate (goat) or anti-human IgA-HRP-conjugate (sheep)	35 ml , ready to use, IgG , coloured red, transparent bottle, red cap IgA , coloured purple, transparent bottle, purple cap
4	SUBSTR TMB	substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	35 ml , ready to use, black bottle, blue cap
5	INCUTRAY	incubation tray with cover	2

Materials required but not provided

- glassware
- pipettes and tips for 15 µl and 1500 µl
- plastic forceps
- rocking platform (vertical)
- filter paper
- collecting devices for infectious material
- distilled or deionized water

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 20 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

Working steps

1. Incubate the test strips with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **15 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **1.5 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **1.5 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **1.5 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.


Test evaluation

Test validity

The *Seraline*[®] Zöliakie-3 IgG / *Seraline*[®] Zöliakie-3 IgA test strips contain 5 consecutive control lines under the strip number:

- Function control (positive reaction with every sample).
- Conjugate controls IgG and IgA (positive reaction with the corresponding conjugate).
Note: The other conjugate control band may develop a weak, non-specific colour.
- The intensity of the negative control must be less than the Cut off control.
- Cut off control (band intensity is used for result evaluation of the diagnostic bands).

The test run is valid if the function control, the corresponding conjugate control and the Cut off control are clearly visible.

Result interpretation

Generally the evaluation of the test has to be done with dried strips. To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet. Alternatively the strips can be evaluated with the *Seraline*[®] scan software.

Judgement of the bands is performed according to the following classification:

Result interpretation	IgG	IgA
negative	Colour intensity of bands < Cut off control	Colour intensity of bands < Cut off control
positive*	Colour intensity of bands ≥ Cut off control	Colour intensity of bands ≥ Cut off control

* A positive detection of anti- *Saccharomyces cerevisiae* antibodies (especially ASCA of isotyp IgA) is an indication for Crohn's disease.

Performance characteristics

Sensitivity

The performance of the *Seraline*[®] Zöliakie-3 IgG / *Seraline*[®] Zöliakie-3 IgA has been compared with another commercially available Line Immunoassay (Test A) by investigation of serum samples from patients with and without clinical suspect of coeliac disease.

Gliadin:

		<i>Seraline</i> [®] Zöliakie-3 IgG	
		positive	negative
Competitor IgG	positive	8	2
	negative	4	22

Agreement: 83%

		<i>Seraline</i> [®] Zöliakie-3 IgA	
		positive	negative
Competitor IgA	positive	6	2
	negative	2	26

Agreement: 89%

Tissue Transglutaminase (tTG):

		<i>Seraline</i> [®] Zöliakie-3 IgG	
		positive	negative
Competitor IgG	positive	5	1
	negative	6	24

Agreement: 81%

		<i>Seraline</i> [®] Zöliakie-3 IgA	
		positive	negative
Competitor IgA	positive	8	0
	negative	2	26

Agreement: 94%

Mannan:

		<i>Seraline</i> [®] Zöliakie-3 IgG	
		positive	negative
Competitor IgG	positive	4	1
	negative	2	29

Agreement: 92%

		<i>Seraline</i> [®] Zöliakie-3 IgA	
		positive	negative
Competitor IgA	positive	3	0
	negative	0	33

Agreement: 100%

Specificity

Results from the investigation of randomly collected blood donor samples (n = 100) using *Seraline*[®] Zöliakie-3 IgG / *Seraline*[®] Zöliakie-3 IgA.
Specificity IgG: 83.0% Specificity IgA: 94.0%

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other Immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only. Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate. Do not use reagents from other manufacturers to complete the kit. Some of the reagents (2, 3, 4) contain small amounts of Kathon (1.0% v/v) as preservative. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use