## INTENDED USE

The Microgen I.M. Kit is intended for *in vitro* detection of heterophile antibodies associated with infectious mononucleosis in human serum or plasma samples. The product is intended for professional use only.

## PRINCIPLE OF THE TEST

The test is based upon the differential absorption test developed by Davidsohn<sup>1</sup> in which horse erythrocytes are agglutinated by heterophile antibodies<sup>2</sup> associated with infectious mononucleosis. Heterophile antibodies to Forssman antigens, which can cause false positive reactions, are absorbed by guinea pig kidney antigen but not by bovine red cell antigen whereas heterophile antibodies to I.M. are absorbed by bovine red cells but not by guinea pig kidney antigen<sup>1,3</sup>. Two samples of the patient's serum or plasma<sup>4</sup> are absorbed with either guinea pig kidney antigen or ox cell antigen on a test slide. Horse red cells are then added and the mixture examined for agglutination. Agglutination in the presence of guinea pig antigen but not with ox cell antigen is indicative of the presence of heterophile antibodies associated with I.M. Other patterns of agglutination are indicative of non-specific reactions and are regarded as negative for I.M.

### CONT

### KIT PRESENTATION

Ag	GPK		
		M121a	2mL

Guinea Pig Kidney Antigen suspended in a solution containing protein stabiliser and 0.5% phenol as preservative. (Yellow dropper)

Ag	OX		
		M121b	2mL

Ox Cell Antigen suspended in a solution containing protein stabiliser and 0.5% phenol as preservative. (White dropper)

REAG	TEST		
		M121c	4mL

Horse Cell Reagent - horse erythrocytes suspended in a stabilising solution containing 0.099% sodium azide as preservative. (Blue dropper)

CONTROL	+		
		M121d	2mL

Positive Control - I.M. positive human serum diluted in a solution containing protein stabiliser and 0.099% sodium azide as preservative. (Black dropper)

Test Slide  
Instructions for Use

## Materials and Equipment required but not provided

- Pipette with disposable tips to deliver 25µL
- Timer
- Mixing sticks (plastic or wood)
- 70% isopropyl alcohol
- Absorbent paper tissues
- Latex gloves, safety glasses and other appropriate protective garments
- Biohazardous waste containers

## WARNINGS AND PRECAUTIONS

### Safety

1. The reagents in this kit are for *in vitro* diagnostic use only.
2. All human blood products should be treated as potentially infectious. The material from which the kit control serum was derived was found negative for antibodies to HIV-1 and HIV-2, HbsAg and HCV when tested. However no known test methods can offer complete assurance that products derived from human blood will not transmit infectious agents.
3. The ox cell antigen is produced from calf blood from animals declared fit for human consumption, and blood components are regarded as having no detectable infectivity for transmissible spongiform encephalopathy agents. Nevertheless, this reagent should be handled as though potentially infectious.
4. The kit components are preserved with either phenol or sodium azide. Although the concentrations are low, take appropriate precautions when handling these reagents. Additionally, sodium azide may react with lead or copper plumbing to form explosive metal azides. Dispose by flushing with a large volume of water to prevent azide build up.
5. Avoid contact of reagents with skin or mucous membranes. If this happens, wash with copious amounts of water.
6. Waste reagents and materials should be decontaminated and discarded according to national guidelines. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30minutes. Liquid waste containing acid must be neutralised before treatment.
7. The test slide should be cleaned thoroughly between samples using 70% isopropanol.

### Procedural

1. The kit should be used in accordance with these Instructions for Use.
2. Use within expiry date marked on the kit.
3. Ensure that all reagents have reached room temperature (18-25°C) before use.
4. Do not mix reagents from different batches.
5. Do not cross-contaminate reagents.
6. The horse cell reagent should always be a dark red colour. Do not use if there is any evidence of cell lysis.
7. Always run the positive control to confirm that the test system is functioning correctly. The positive control has been diluted so that it will only give a weak reaction. Should the kit be incorrectly stored or used beyond its expiry date, the control may give a negative reaction. Under such circumstances, the reagents should be discarded and a fresh kit used.

## STORAGE AND SHELF LIFE

### Kit components:

All unopened components can be used until the date printed on the outer carton label provided they are stored at 2-8°C. Once opened, reagents can be used repeatedly during the stated shelf life of the kit provided that they are stored at 2-8°C and that reasonable precautions are taken to avoid high levels of microbial contamination.

### Specimens:

Non-haemolysed serum or plasma can be used. Inactivated serum is also suitable. Separated samples should be stored at 2-8°C and tested within 2-3 days if possible. If longer-term storage is necessary, specimens should be frozen at -20°C or below. Specimens that have been frozen and then thawed should be thoroughly mixed before testing.

## TEST PROCEDURE

1. Ensure that all reagents are at room temperature (18-25°C) and mixed thoroughly.
2. Carry out the following procedure on the white glass test slide provided. Always run the positive control to confirm the test system is functioning correctly.

Oval A	Oval B
Dispense 25µL patient's sample or positive control.	Dispense 25µL patient's sample or positive control.
Add 1 drop guinea pig kidney antigen (Yellow dropper)	Add 1 drop ox cell antigen (White dropper)
Mix and allow to stand for 30 seconds	Mix and allow to stand for 30 seconds
Add 1 drop horse cell reagent (Blue dropper)	Add 1 drop horse cell reagent (Blue dropper)
Mix and rock slide gently for 1 minute	Mix and rock slide gently for 1 minute

3. Examine for agglutination at **1 minute**.

## INTERPRETATION

1. A **positive** result is indicated by distinct agglutination in test oval A (guinea pig antigen absorption) with no significant agglutination in test oval B (ox cell absorption). This is a true **differential absorption** as demonstrated by the positive control.
2. The result is **negative** when there is no discernible agglutination in either test oval A or B.
3. Any sample which does not present a clear result as defined above should be diluted 1:5 in 0.85% sodium chloride solution and the test repeated. Criteria 1 and 2 above should be applied to the repeat test result.
4. The positive control (tested undiluted) should give a weak positive reaction in Oval A with no appreciable agglutination in Oval B. If no agglutination is seen in Oval A, the test is invalid as the reagents have probably deteriorated.
5. The test is invalid if any result other than those defined in 1-4 above is obtained.

## LIMITATIONS OF USE

1. Results from the Microgen I.M. Kit should always be interpreted in the context of all available clinical and laboratory data.
2. If the blood film picture does not agree with the Microgen I.M. Kit results, the test should be repeated on a fresh sample drawn one week later. Testing for evidence of EB virus may also aid diagnosis.
3. There are some non-EBV associated pathological conditions which can produce false positive results in haemagglutination tests of this type. These include some other viral infections and lymphomatous conditions.
4. The use of EDTA plasmas, rather than serum, can also be associated with a higher incidence of potentially misleading results. However, careful observance of the interpretation instructions above will minimise erroneous conclusions.

## PERFORMANCE CHARACTERISTICS

Microgen I.M. Kit has been independently evaluated in two leading microbiology laboratories in the United Kingdom.

### Laboratory 1

As part of an evaluation of 14 commercial test kits, 232 plasma and sera were tested by the Microgen I.M. Kit in parallel with a reference system comprising an immunofluorescent assay for IgM antibodies to EBV viral capsid antigen, backed up by haematology profile.

Table 1: Comparison of Microgen I.M. Kit with reference IFA

	Microgen I.M. Kit		Total
	Positive	Negative	
Reference Methods	73	0	73
	6	153	159
Total	79	153	232

Sensitivity	100%
Specificity	96.2%
Accuracy	97.4%

Six false positive results were seen:

- a) 2 specimens were from patients who presented as suspected I.M. but reference methods (IFA, blood film and blood counts) were negative.
- b) 1 patient presented with lymphadenopathy
- c) 1 patient presented with chronic fatigue
- d) 2 samples gave very weak agglutination. No clinical symptoms were available although the samples had been sent by the G.P. for I.M. screening.

### Laboratory 2

As part of an evaluation of 6 commercially available test kits, 40 plasma specimens received by the laboratory for routine Full Blood Count and I.M. screening were tested by the Microgen I.M. Kit. A well-established and respected commercial assay was used as the reference procedure in conjunction with microscopic examination of blood films for the presence of atypical lymphocytes, a characteristic of infectious mononucleosis.

Table 2: Comparison of Microgen I.M. Kit with reference assay

	Microgen I.M. Kit		Total
	Positive	Negative	
Reference Assay	24	0	24
	0	16	16
Total	24	16	40

Sensitivity	100%
Specificity	100%
Accuracy	100%

In this trial no discrepant results were obtained.

## REPRODUCIBILITY

**Intra-batch reproducibility** – This has been evaluated by testing one batch of kits on 3 different occasions using 2 operators. On each occasion the sensitivity of the kit was tested with serial dilutions of standard reference and kit control sera. No difference in end-point titres were seen between different operators/different test occasions

**Inter-batch reproducibility** – Comparison of the QC release data for 30 batches of I.M. Kit has shown that there is only marginal variation of sensitivity between batches (maximum 1 doubling dilution of reference positive serum). Inter-batch reproducibility of specificity was confirmed by testing 3 batches of kits with a panel of 15 randomly selected plasma specimens. Identical results were obtained with all batches.

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**BIBLIOGRAPHY / BIBLIOGRAFÍA / BIBLIOGRAPHIE**

1. Davisohn, I., Stern, K. and Kashwagi, C. "The Differential Test for Infectious Mononucleosis." Amer. J. Clin. Path., 21 1101 (1951)
2. Paul, J.R. and Bunnell, W.W. "The Presence of Heterophile Antibody in Infectious Mononucleosis." Amer. J. Med., 183 909 (1932)
3. Lee, C.L., Davisohn, I. and Slaby, R. "Horse Agglutinins in Infectious Mononucleosis." Amer. J. Clin. Path., 49 3 (1968)
4. Davidson, R.J.L. and Main, S.R. "Use of Plasma instead of Sera in Laboratory Test for Infectious Mononucleosis". J. Clin. Path., 24 259 (1971)

**Microgen Bioproducts Ltd**  
**1 Admiralty Way**  
**Camberley**  
**Surrey**  
**GU15 3DT**

**WF 0960/2003-03**