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


REFERENCES

1. Marshall BJ, Warren JR. Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. Lancet 1984;i: 1311-1314.
2. Versalovic J. *Helicobacter pylori*. Pathology and diagnostic strategies. Am J Clin Pathol 2003;119:403-12.
3. Dunn BE, Cohen H, Blaser MJ. *Helicobacter pylori*. Clin. Microbiol. Rev. 1997;10:720-41.
4. Gregson DB, Simor AE. *Helicobacter pylori* and inflammatory disease of the stomach and duodenum. Stomach and Bowel, 1991;Sep:29-30, 35-7.
5. Parsonnet J, Friedman GD, Vandersteen DP et al. 1991. *Helicobacter pylori* infection and the risk of gastric carcinoma. N Engl J Med 1991;325:1127-31.
6. Sanders MK, Peura DA. *Helicobacter pylori*-associated diseases. Curr Gastro Reports. 2002;4:448-54.
7. Uemura N, Okamoto S, Yamamoto S et al. *Helicobacter pylori* infection and the development of gastric cancer. N Engl J Med. 2001;345:784-9.
8. Graham DY. I. *Helicobacter pylori*: Its epidemiology and its role in duodenal ulcer disease. J Gastroenterol Hepatol 1991;6:105-13.
9. Malaty HM, Graham DY, Klein PD et al. Transmission of *Helicobacter pylori* infection. Studies in families of healthy individuals. Scand J Gastroenterol 1991;26:927-32.
10. Marshall BJ, Armstrong JA, McGeachie DB, Glancy RJ. Attempt to fulfill Koch's postulates for pyloric campylobacter. Med J Austral 1985;142:436-439.
11. Veralovic J, Fox JG. *Helicobacter*. In: Murray PR et al eds. Manual of clinical microbiology, 8th ed. Washington DC:ASM Press, 2003:915-28.

12. Alpert LC, Graham DY, Evans DJ Jr et al. Diagnostic possibilities for *Campylobacter pylori* infection. Eur J Gastroenterol Hepatol 1989;1:17-26.
13. Barthel JS, Everett ED. Diagnosis of *Campylobacter pylori* infections: The "gold standard" and the alternatives. Rev Infect Dis 1990;12:S107-14.
14. Nichols L, Sughayer M, DeGirolami PC et al. Evaluation of diagnostic methods for *Helicobacter pylori* gastritis. Am J Clin Pathol 1991;95:769-73.
15. Vaira D, Vakil N. Blood, urine, stool, breath, money and *Helicobacter pylori*. Gut 2001;48:287-9.
16. Graham DY, Klein PD, Evans DJ Jr et al. *Campylobacter pylori* detected noninvasively by the 13C-urea breath test. Lancet 1987;i:1174-7.
17. Graham DY, Malaty HM, Evans DG et al. *Epidemiology of Helicobacter pylori* in an asymptomatic population in the United States. Effect of age, race, and socioeconomic status. Gastroenterology 1991;100:1495-1501.
18. Talley NJ, Newell DG, Ormand JE et al. Serodiagnosis of *Helicobacter pylori* Comparison of enzyme-linked immunosorbent assays. J Clin Microbiol 1991;29:1635-9.
19. Vaira D, Malfertheiner P, Mégraud F et al. Diagnosis of *Helicobacter pylori* infection with a new non-invasive antigen-based assay. Lancet 1999;354:30-3.
20. Vaira D, Malfertheiner P, Mégraud F et al. Non invasive antigen-based assay for assessing *Helicobacter pylori* eradication: a European Multicentre Study. Am J Gastroenterol 2002;95:925-9.
21. Montiero L, de Mascarel A, Sarraqueta A. et al. Diagnosis of *Helicobacter pylori* Infection: Noninvasive Methods Compared to Invasive Methods and Evaluation of Two New Tests. Am J Gastroenterol 2001;96:353-8.
22. Malfertheiner P, Megraud F, Morain O et al. Current concepts in the management of *Helicobacter pylori* infection – The Maastricht 2-2000 Consensus Report. Aliment Pharmacol Ther 2002;16:167-80.
23. Murray DM. Clinical relevance of infection by *Helicobacter pylori*. Clin Microbiol Newlett 1993;15:33-7.
24. Loffeld R.J.L.F., Stobberingh E, Van Spreeuwel JP et al. The prevalence of anti-*Helicobacter (Campylobacter) pylori* antibodies in patients and healthy blood donors. J Med Microbiol 1991;32:105-9.
25. Morris A, Nicholson G, Lloyd G et al. Seroepidemiology of *Campylobacter pyloridis*. N. Z. Med. J. 1986;99:657-9.
26. Fiedorek SC, Malaty HM, Evans D Letal. Factors influencing the epidemiology of *Helicobacter pylori* infection in children. Pediatrics 1991;88:578-582.

SN11269

REV. 05/12

| | |
|---|--|
|  Manufactured For | Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: (513) 271-3700 Orders/Customer Service: (800) 543-1980 Technical Support Center: (800) 343-3858 Information Fax:(513) 272-5432 Ordering Fax:(513) 271-0124 |
| | Meridian Bioscience Europe Via dell'Industria, 7 20020 Villa Cortese (MI) Italy Tel: +39 0331 433636 Fax: +39 0331 433616 e-mail: info@mdeur.com |

Meridian Bioscience Europe s.a./ n.v.
 Rue de l'Industrie 7
 1400 Nivelles
 Belgium
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 e-mail: info@mdeur.be

Meridian Bioscience Europe b.v.
 Halderheiweg 6
 5282 SN Boxtel
 The Netherlands
 Tel: +31 (411) 621166
 Fax: +31 (411) 624841
 e-mail: meridian.info@planet.nl

Meridian Bioscience Europe France
 34, rue de Ponthieu
 75008 Paris
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 e-mail: info@meridianbioscience.fr

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:
Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Erläuterung der graphischen symbole)

| | | | |
|---|--|---|---|
|  | Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis | CONTROL + | Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle |
| LOT | Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung | CONTROL - | Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle |
| IVD | In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum | EC REP | Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft |
|  | This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG. | SMP PREP DIL SPE | Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet |
| REF | Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer |  | Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren |
|  | Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten | BUF RXN | Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer |
|  | Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller |  | For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung |
|  | Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen | SOLN STOP | Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung |
|  | Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung | CONJ ENZ | Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat |
| SN | Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer | CONTROL | Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest |
| TEST | Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät | REAG | Reagent / Reagente / Réactifs / Reactivos / Reagenzien |
|  | Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum | BUF WASH | Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer |
| BUF | Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer |  | Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise |
| CONJ | Conjugate / Coniugato / Conjugué / Conjugado / Konjugat | DIL SPE | Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillon / Diluyente de muestra / Probenverdünnungspuffer |
| SUBS | Substrate / Substrato / Substrat / Substrato / Substrat | BUF WASH 20X | Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat |

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.