## REACTIVIDAD CRUZADA

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### REFERENCES

- Bartlett JG. *Clostridium difficile*: clinical considerations. Rev Infect Dis. 1990;12 S243-S251.
- Bartlett JG, TW Chang, M Gurwith, SL Gorbach and AB Onderdonk. Antibiotic associated pseudomembranous colitis due to toxin producing clostridia. N Engl J Med 1978;298: 531-534.
- George WL, VL Sutter, EJC Goldstein, SL Ludwig and SM Finegold. Aetiology of anti-microbial agent associated colitis. Lancet. 1978;1: 802-803.
- Rolle RD. Asymptomatic intestinal colonization by *Clostridium difficile*. In: RD Rolfe and SM Finegold (eds.). *Clostridium difficile*: its role in intestinal disease. Academic Press, Inc. New York. 1988.
- McFarland LV, ME Mulligan, RYY Kwok and WE Stamm. Nosocomial acquisition of *Clostridium difficile* infection. New Engl J Med. 1989;320:204-210.
- Johnson S, A Adelman, CR Clabots, LR Peterson, DN Gerding. Recurrences of *Clostridium difficile* diarrhea not caused by the original infecting organism. J Infect Dis. 1989;159 (2):340-343.
- Sullivan NM, S Pellett, TD Wilkins. Purification and characterization of Toxins A and B of *Clostridium difficile*. Infect Immun. 1982;35(3): 1032-1040.
- Lyerly DM, CJ Phelps, J Toth and TD Wilkins. Characterization of toxins A and B of *Clostridium difficile* with monoclonal antibodies. Infect Immun. 1986;54(1):70-76.
- Lyerly DM, CJ Phelps, and TM Wilkins. Monoclonal and specific polyclonal antibodies for immunoassay of *Clostridium difficile* Toxin A. J Clin Microbiol. 1985;21(1): 12-14.
- Chang TW, M Lauerhmann and JG Bartlett. Cytotoxicity assay in antibiotic-associated colitis. J Infect Dis. 1980;140:765-770.
- Lyerly DM, KE Saum, DK McDonald, and TD Wilkins. Effects of *Clostridium difficile* toxins given intragastrically to animals. Infect Immun. 1985;47:349-352.
- Meador J, and RK Tweten. Purification and characterization of toxin B from *Clostridium difficile*. Infect Immun. 1988;56(7):1708-1714.
- Data on file, Meridian Bioscience, Inc.
- Cooperstock, M. *Clostridium difficile* in infants and children. In: RD Rolfe and SM Finegold (eds.). *Clostridium difficile*: its role in intestinal disease. Academic Press, Inc. New York. 1988.
- Lyerly DM and TD Wilkins. Purification and properties of toxins A and B of *Clostridium difficile*. In: RD Rolfe and SM Finegold (eds.). *Clostridium difficile*: its role in intestinal disease. Academic Press, Inc. New York. 1988.
- Willey SH and JG Bartlett. Culture for *Clostridium difficile* in stools containing a cytotoxin neutralized by *Clostridium sordellii* antitoxin. J Clin Micro. 1979;10: 860-864.

SN11161

REV. 01/15

 Manufactured By	Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
 Authorized Representative	Meridian Bioscience Europe S. r. L Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.  
 2 Avenue du Japon - 1420 Braine l'Alleud  
 BELGIUM  
 Tel: +32 (0) 67 89 59 59  
 Fax: +32 (0) 67 89 59 58  
 Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France  
 34 rue de Ponthieu - 75008 Paris  
 FRANCE  
 Tel: +33 (0) 1 42 56 04 40  
 Fax: +33 (0) 9 70 06 62 10  
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.  
 Postbus 301 - 5460 AH Veghel  
 NETHERLANDS  
 Tel: +31 (0) 411 62 11 66  
 Fax: +31 (0) 411 62 48 41  
 Email: info.bn@meridianbioscience.eu

### INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
<b>CE</b>	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostico in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In-Vitro-Diagnostika	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenverarbeitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero di catalogo / Référence ou catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contents sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contenu suffisant pour <n> tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Límites de temperatura / Temperaturbegrenzung	<b>CONJ ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Coniugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
<b>BUF</b>	Buffer / Soluzione tamponne / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Coniugato / Coniugé / Coniugado / Konjugat	<b>DIL SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF WASH 20X</b>	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		<b>DET REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagent

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.