



QuantideX[®]
qPCR DNA QC Assay

MEASURE THE DNA THAT MATTERS.

ASSAY FEATURES

- ▶ MEASURES THE DNA THAT MATTERS: PCR-AMPLIFIABLE DNA
- ▶ PROVIDES ACTIONABLE GUIDANCE THAT RESCUES LOW-QUALITY DNA AND HELPS ASSURE THE ACCURACY OF NGS CALLS
- ▶ A 2-IN-1 ASSAY: DISTINGUISHES LOW-QUALITY, UNAMPLIFIABLE DNA FROM DNA WITH SAMPLE-DERIVED PCR INHIBITION
- ▶ SIMPLE, HIGH-THROUGHPUT WORKFLOW USING A CONVENIENT MULTIPLEXED DESIGN

EXTENSIVELY VALIDATED IN PEER-REVIEWED PUBLICATIONS

1866 RESIDUAL CLINICAL FFPE SPECIMENS FROM 13 COHORTS




Figure 1 : The QuantideX® qPCR DNA QC Assay has been assessed with >2000 FFPE and FNA tumor biopsies across more than a dozen different cohorts. Each sample is represented by its corresponding QFI™ Score.

QuantideX[®]

qPCR DNA QC Assay

GUIDES DNA INPUT TO IMPROVE THE ACCURACY OF NGS VARIANT CALLS AND RESCUES FALSE POSITIVE AND NEGATIVE CALLS




Figure 2 : The specificity of NGS mutation detection can be dramatically improved by increasing the input of amplifiable DNA from a low-quality FFPE tumor biopsy (QFI Score= 2.6%). Although the expected BRAF V600E mutation is detected at all inputs, the number of false positives is significantly increased at lower copy number inputs. The QuantideX[®] qPCR DNA QC Assay informs an input of 200 amplifiable copies, which is associated with detection of a single variant, the expected BRAF V600E mutation. The mass inputs at the top of each bar were measured using spectrophotometry.

IDENTIFIES PCR INHIBITORS PRESENT IN THE SAMPLE




Figure 3 : Data from 15 melanoma FFPE DNA samples isolated with three different isolation methods (IM1, IM2, IM3). A high Cq is an indicator of inhibition in a sample.

PROVIDES ACTIONABLE GUIDANCE FOR DNA INPUTS INTO DOWNSTREAM NGS ENRICHMENT STEPS

SAMPLE NAME	DNA INPUT (ng)	INHIBITION	AMPLIFIABLE COPY NUMBER PER μ L	QFI [™] SCORE
SAMPLE A	10	PASS	300	10.0
SAMPLE B	10	PASS	856	28.5
SAMPLE C	15	FAIL	FAIL	FAIL
SAMPLE D	5	AT RISK	56	3.7
SAMPLE E	20	PASS	8	0.1

Figure 4 : Actionable results generated from the QuantideX[®] qPCR DNA QC Assay. The results identify PCR inhibition, reveal the DNA quality score (QFI[™] Score), and prescribe volume inputs based on sample-specific DNA copy number for theoretical detection of a 5% variant.

KIT ORDERING INFORMATION

QUANTIDEX® qPCR DNA QC ASSAY
[P/N 49539] 100 Reactions.

*Research Use Only. Not for use in diagnostic procedures.

REFERENCES

1. Functional DNA quantification guides accurate next-generation sequencing mutation detection in formalin-fixed, paraffin-embedded tumor biopsies. Sah et al. *Genome Med* 2013;5(8):77
2. Evaluation of an integrated clinical workflow for targeted next-generation sequencing of low-quality tumor DNA using a 51-gene enrichment panel. Choudhary et al. *BMC Med Genomics* 2014;(7):62
3. Variation in pre-PCR processing of FFPE samples leads to discrepancies in BRAF and EGFR mutation detection: A diagnostic RING trial. Kapp et al. *J Clin Pathol* 2015;(68):111-118

Vertrieb durch:

VIROTECH Diagnostics GmbH
Löwenplatz 5
D-65428 Rüsselsheim
☎ +49 (0)6142 69090



Asuragen®
Oncology

Asuragen, Inc. | 2150 Woodward St., Suite 100 | Austin, TX 78744 USA

asuragen.com

2500-0495