

# LINE Software LabImage® Line Analysis (LA)

**Interpretation Software for the  
Virotech LINE Product Range**

ROUTINE SOLUTION AT THE  
HIGHEST AUTOMATION LEVEL




**EASY ---- QUICK ---- RELIABLE**




# A FEW CLICKS TO EVALUATE STRIP BASED ASSAYS





## I



## II



## III



### INTUITIVE AND EASY TO USE:

- ▼ Intuitive guidance through the software
- ▼ User-friendly interface
- ▼ Fast creation of a working list or
- ▼ Full connection to LIMS


### SEVERAL WAYS TO SCAN:

- ▼ Flatbed scanner
- ▼ LabImager TR reader w/wo integrated PC
- ▼ Excellent stability of trays in the reader

### FLEXIBILITY AND INDIVIDUALITY FOR YOUR LABORATORY:

- ▼ Adaptable user interface to create your working list
- ▼ Manual entry possible
- ▼ Designing laboratory-specific reports


# HIGHEST RELIABILITY FOR YOUR ANALYSIS



Strip Detection Result

- ▼ Innovative technology for accurate detection of the strips
- ▼ Defined detection window for reliable detection of the bands
- ▼ Immediate control on scanning: correctly detected strips are marked in green

Detection Notes  
**Number of detected strips: 13**



## EFFICIENT EVALUATION:

- ▼ Single view of each strip for a better monitoring of the results
- ▼ Manual corrections are possible
- ▼ Proceeded corrections are clearly marked – retraceability ensured
- ▼ Reduction of background possible

## FULLY AUTOMATED STRIP ANALYSIS

- ▼ Automatic identification of antibody classes
- ▼ Joined analysis of IgG- and IgM-results (e.g. for Borrelia)
- ▼ **Several antibody classes in one run** (e.g. IgM and IgG)
- ▼ **Various parameters in one run** (e.g. Borrelia and Treponema)
- ▼ Assistance in every aspect


Name	Surname	IgG Evaluation	IgM Evaluation	Interpretation
Adrian	Armano	borderline	positive	positive / suspicious for early Lyme Borreliosis
Bertha	Bell	negative	borderline	borderline / potential indication for Early Lyme Borreliosis, control recommended
Charles	Chaplin	negative	positive	positive / suspicious for early Lyme Borreliosis
David	Hassle	negative	negative	negative / no indication of Lyme Borreliosis
Eva	Ericson	negative	negative	negative / no indication of Lyme Borreliosis
Felix	Fallon	positive	positive	positive / suspicious for deseminated Lyme Borreliosis
Geraldine	Gardin	positive	positive	positive / suspicious for deseminated Lyme Borreliosis
Henry	Holder	negative	negative	negative / no indication of Lyme Borreliosis
Irvine	Ildig	negative	negative	negative / no indication of Lyme Borreliosis
Joseph	Jordan	borderline	positive	positive / suspicious for early Lyme Borreliosis

## SECURE ARCHIVING OF RESULTS

- ▼ Save the evaluation in a database
- ▼ Possibility to save strips and original tray images
- ▼ Enduring batch documentation
- ▼ Patient history recallable
- ▼ Patient-specific analysis exportable

**The Virotech SetCards provide maximum safety.  
Software and Virotech SetCards are CE-marked  
according to the current IVDD (98/79/EC)!**


# General Overview




**Search:**  
The search function enables the representation of test- and / or patient-related information.



**NEW:**  
A project can be created. According to the project name and its location the following options are available:  
1. Scan  
2. Create work list  
3. Import from LIMS



**Managing:**  
The administration provides access to all saved projects. For sorting folders or unique analysis names can be created.



**After scanning in the detection**  
1. Correction of patient data  
2. Import of parameters and controls



Problemnr.	Typ	Testtyp	Sp-X	Vorname	Nachname	Geburtsdag	SpZ Bereich
1	PAZ	Brennstoff	SP5	Andreas	Andreas	02.02.2000	regional
2	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
3	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
4	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
5	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
6	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
7	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
8	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
9	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
10	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
11	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
12	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
13	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
14	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
15	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
16	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
17	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
18	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
19	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
20	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
21	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
22	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
23	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional
24	PAZ	Brennstoff	SP5	Christoph	Christoph	02.02.2000	regional



**Possible actions:**  
1. Set bands invalid  
2. Detailed view  
3. Edit test- / patient-related information  
4. Export to LIMS  
5. Create reports

## EASY ---- QUICK ---- RELIABLE

LabImage® LA is a sophisticated software for efficient strip recognition and interpretation. The LabImage® LA has an intuitive interface to prevent unnecessary clicks. The complete automation of the strip analysis offers maximum standardisation, which significantly contributes to increase efficiency in medical diagnostic laboratories.

## ORDER INFORMATION

### Order number:

F/LabImage1  
F/LabImage2  
F/LabImage3

### Description:

LabImage® Package 1  
LabImage® Package 2  
LabImage® Package 3

Products	Package 1	Package 2	Package 3
LabImage® Reader with PC	X		
LabImage® Reader without PC		X	
Flatbed Scanner (Plustek)			X
LabImage® Software (CD)	X	X	X
Virotech SetCards (CD)	X	X	X
Manuals	X	X	X
License Dongle		X	X

## VIROTECH LINE IMMUNOBLOT RANGE

### LINES:

▼ B.pertussis + CatACT LINE	IgG/IgA
▼ Borrelia LINE	IgG/IgM
▼ Borrelia in vivo LINE plus	IgG
▼ Borrelia Europe LINE	IgG/IgM
▼ Borrelia Europe plus TpN17 LINE	IgG
▼ EBV LINE	IgG/IgM
▼ Helicobacter pylori LINE	IgG/IgA
▼ HSV LINE	IgG
▼ Mycoplasma pneumonia LINE	IgG/IgM/IgA
▼ Treponema pallidum LINE	IgG/IgM
▼ Yersinia enterocolitica LINE	IgG/IgA

### BENEFITS OF OUR VIROTECH LINE TEST SYSTEM:

- ▼ Solid strips stabilised on plastic foil
- ▼ Assorted strips in a booklet
- ▼ Clear strip numbering
- ▼ Quality control on each strip
- ▼ Identical incubation times (30'/30'/10')
- ▼ Test results after approx. 100 minutes
- ▼ General components (EcoSystem)
- ▼ Cut-off control included in the kit, also available separately
- ▼ Easy reading & interpretation
- ▼ Positive and negative controls separately available

## YOUR CONTACT AT VIROTECH DIAGNOSTICS GMBH

### **Andreas Beczkowski**

Software Specialist

Sales Operation

Mobile: +49 (0)1761 0213465

[andreas.beczowski@virotechdiagnostics.com](mailto:andreas.beczowski@virotechdiagnostics.com)

### **Michael Becker**

Sales Director Export

Tel: +49 (0)6142 690931

Fax: +49 (0)6142 966613

Mobile: +49 (0)1726744920

[michael.becker@virotechdiagnostics.com](mailto:michael.becker@virotechdiagnostics.com)

**Kapelan Bio-Imaging GmbH** is the **Manufacturer** of LabImage® software LA,  
LabImager TR (LR-1000) and LabImager TR with PC (LR-1010)



B\_LabImage\_20170125-02e



VIROTECH Diagnostics GmbH Löwenplatz 5 D-65428 Rüsselsheim

+49 (0)6142 69090 +49 (0)6142 690919

info@virotechdiagnostics.com www.virotechdiagnostics.com

