EurobioPlex

*Bartonella spp/Bartonella henselae/
et Bartonella quintana*

PCR EN TEMPS REEL

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Version 1.03 from February 2nd 2016

Storage conditions:

Keep all reagents at -20°C until use and after first use

Instructions for use

INTENDED USE

The Eurobioplex Bartonella EBX-013 test uses real time polymerase chain reaction (PCR) amplification and is designed for the qualitative detection of Bartonella spp. (all Bartonella), Bartonella henselae and Bartonella quintana from a DNA nucleic acid extract. This kit also contains a control of inhibition of PCR. The test allows early diagnosis of presumption of Bartonella infection in humans, from serum, blood or any biological samples of patients, or complement a proven or indeterminate serological diagnosis by standard techniques (serology, histology and culture).

INTRODUCTION

Bartonella are bacteria that infect many species of mammals, including humans. Chronic infection of erythrocytes is a specific adaptation of Bartonella to mammals. They can be also found in the skin, the bone tissue and endothelial cells. The main species of Bartonella described in the literature have been isolated in the 1990s. About half of them are implicated as agents of zoonoses. The major reservoir is then domestic and wild animals living close to humans. For example in France, it is mainly cats. Human diseases caused by Bartonella spp. are called the Bartonellosis. The two main Bartonella species are Bartonella henselae responsible for “cat scratch disease”, and B Quintana responsible for “trench fever” and endocarditis. Other less common species of Bartonella are likely to cause other human pathologies: Bartonella bacilliformis, B vinsonii berkhoffi, B clarridgeiae, B tamiae, B rochalimae, B elizabethae, B koehlerae, B graham and B alsatica. Endocarditis are the most serious complication.

Depending on the species Bartonella have various modes of transmission. Humans can be infected with the animal directly by scratch, bite, or simple caresses if the skin of the hands is scratched. The infection can also be transmitted indirectly through a blood-sucking arthropod vector. Each species of Bartonella usually has a single natural reservoir host which is generally asymptomatic due to a long-term bacteremia. The accidental hosts however develop symptoms generally more severe. It is difficult to propose a treatment for bacteria causing so many variable clinical manifestations. Moreover, the evolution of the disease varies according to the bacteria and immune system of the host. It can be acute, recurrent or even chronic. The treatment has to be adapted to each clinical situation. There is so far no vaccine prophylaxis.

Cat scratch disease (CSD)

The pathogen agent is Bartonella henselae. The CSD is characterized by infectious conjunctivitis accompanied by a strong lymph node reaction. It was described for the first time in 1889 by H.Parinaud. The prevalence of the CSD is 6.6 cases per 100,000 population. The peak of incidence is reported during the cold period of the year, between September and April when pets live more in contact with their master at home.

The typical symptoms are lymphadenopathy, fever (38-41°C), fatigue, headaches, splenomegaly and cutaneous manifestations such as transitional maculo-papular rash. The disease vectors are fleas with a contamination by the way of a scratch or a bite of an infected cat or through the bite of an infected flea. It is necessary to adopt a number of preventive measures for exposure to animals, especially domestic pets.

The prognosis is good in immunocompetent patients, with spontaneous resolution in 90% of cases. On the other hand, antibiotic therapy is often necessary for children and immunocompromised people, as well as in cases of neurological and cardiac forms.

Trench fever

The pathogen agent is *Bartonella Quintana*. It is an intracellular and ubiquitous germ transmitted by body lice causing periodic fever with headaches and severe pain of the lower limbs. Fever trenches is the clinical principal form of this infection that has been essentially described during the First World War causing then epidemics with more than a million cases estimated in Europe. It reappeared during the Second World War and finally disappeared. Currently it is a re-emerging disease which primarily concerns populations exposed to precarious hygiene and the presence of body lice. It notably affects homeless persons in the industrialized countries with a high prevalence (16 to 30%).

Diagnosis

Diagnosis of these infections relies on serology, histology, blood culture and the research of bacterial DNA from different samples. *Bartonella* are small gram-negative bacilli, aerobic, polymorphic, non-acido-alcohol-resistant, often immobile, from 1 to 1,2µm long and 0,5µm in diameter. They are slow-growing, as difficult to cultivate to isolate. However, these bacteria can be grown on media with fresh blood, humid atmosphere enriched in CO₂, incubated at 35°C, but growth of bacteria on these environments remains particularly tedious and costly. Diagnosis is based first on the clinical description and then on serology (indirect immunofluorescence or immuno-assay enzyme with titration of IgG or IgM antibodies). The use of the techniques of molecular biology (PCR and sequencing) allows a rapid diagnosis from tissues (lymph node, skin, heart valve,...). It also enables the diagnosis of species. The diagnosis by qPCR brings sensitivity, reproducibility and early detection of the presence of the bacteria.

PRINCIPLE OF THE DETECTION

The Eurobioplex EBX-013 test is a test of amplification of nucleic acid of *Bartonella* spp., *Bartonella henselae* and *Bartonella quintana* by real-time PCR as well as a control of PCR inhibition. This control ensures that a negative result may not be due to the presence of PCR inhibitors at high quantity.

The test is performed from the DNA extracted from the sample. *Bartonella* spp, *henselae* and *quintana* are respectively detected using a FAM, HEX, and TEXAS RED labeled probe. The specific probe of the PCR inhibition control is labelled by the fluorochrome CY5. All emit a specific fluorescence following its hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence relates the accumulation of specific amplification products.

This amplification system has been validated on:

- a plasmid containing specific inserts of *Bartonella* spp., *Bartonella henselae* and *Bartonella quintana* (positive control: plasmid CP-BsppHQ),
- a non-coding random sequence inserted into plasmid PUC57 for the control of PCR inhibition (plasmid CI-PCR),
- clinical samples of known contamination (see clinical study).

DESCRIPTION AND CONTENT OF THE KIT

Bartonella spp/henselae/quintana real-time PCR kit is ready to use for the specific detection of bacteria belonging to this genus and species. The kit contains reagents and enzyme necessary for the amplification of Bartonella spp, Bartonella henselae and Bartonella quintana DNA, and of the plasmid control of PCR inhibition. Fluorescence is emitted and measured individually by an optical system during PCR. The detection of the amplified fragment is carried out by a fluorimeter using the channels listed in the table 1 below:

Table 1 :

Target	Fluorophore	Excitation	Emission
Bartonella spp	FAM	495 nm	515 nm
Bartonella henselae	HEX	535 nm	555 nm
Bartonella quintana	TexasRed	585 nm	605 nm
Control of PCR inhibition	Cy5	647 nm	667 nm

Alternatives:

- Channel **FAM** (ABI Systems, SmartCycler II, Chromo 4/CFX96, Systems Mx), Channel 530 (LC 480), Channel Green (RotorGene)
- Channel **HEX** (Chromo 4, Systems Mx), Channel VIC (ABI Systems, Biorad CFX96), Channel Alexa532 (SmartCycler II), Channel 560 (LC 480), Channel Yellow (RotorGene),
- Channel **Texas Red** : LC Red 610 (LC480), Channel Orange (Rotor Gene)
- Channel **Cy5** (ABI Systems, Chromo 4/CFX96, Mx Systems), Channel 670 (LC480), Channel Alexa647 (SmartCyclerII), Channel Red (RotorGene)

Table 2 :

Components of the kit	96 tests	48 tests	24 tests	Reconstitution
Master Mix	4x530 µl	2x530 µl	530 µl	Ready to use
Water = negative control (CN)	1 ml	1 ml	1 ml	Ready to use
Positive control (CP-BsppHQ)	50 µl	50 µl	50 µl	Ready to use
Control of PCR inhibition	120 µl	60 µl	30 µl	Ready to use

Required material not provided:

- ◇ Biological Hood
- ◇ qPCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for qPCR
- ◇ Micropipettes
- ◇ Filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves

STORAGE

All reagents must be stored at -20°C. Storage at +4°C is not recommended.
All reagents can be used until the expiration date indicated on the kit.
Many freezing/defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

CAUTIONS AND NOTES

Read carefully instructions before starting.

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ Many freezing/defrosting cycles of the reagents must be avoided, and could lead to decrease in sensitivity.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ Use of ice or cooling block is recommended in case of delays due for example to a significant number of samples to treat or high temperatures.
- ◇ Define three working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures. The use of gloves is mandatory. Appropriate methods of preparation of DNA should be used.
- ◇ Always use filtered tips for micropipettes.
- ◇ Use specific lab coat and gloves in each working area.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

COLLECTION OF SAMPLES, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes.
- ◇ The samples should be extracted immediately or frozen from -20°C to -80°C.
- ◇ Transport of clinical samples must obey local regulations for this type of infectious agents.

PROCEDURE

I- Extraction of the DNA

DNA extraction kits are available from several manufacturers. You can use your own homemade technique or a suitable commercial system by referring to the manufacturer's instructions.

For liquid samples, InnuPREP Blood DNA Mini kit (845-KS-1040050, from Eurobio) is recommended following the protocol for extraction of body fluids in a final volume of elution of 50 μ L.

The control of PCR inhibition (CI-PCR; channel Cy5) added to the PCR reaction ensures that a negative result may not be due to the presence of PCR inhibitors at high quantity. It can also be added before the extraction and makes sure that a negative result may not be due to a problem of extraction or the presence of PCR inhibitors at high quantity. CI-PCR is also available from Eurobio (Ref EurobioPlex EBX-002). In the 2 cases we recommend the addition of 1 microlitre per extraction (in a final volume of elution of 50 μ L) or by PCR reaction.

II- qPCR

General comments:

- a) Homogenize the reaction mix before starting.
- b) Positive controls (BspHQ and control of PCR) contain a very high concentration of DNA. Manipulations must be performed carefully to avoid contamination.
- c) Dilutions may be kept maximum 24 hours at + 4°C.

It is necessary to test: a positive control CP-BspHQ (Bartonella spp/ henselae/ quintana), as well as a negative control (water supplied) and the control of PCR inhibition (provided control = CI-PCR). No standard range is performed.

Example of qualitative qPCR plate layout:


	1	2	3	4	5	6	7	8	9	10	11	12
A	CP-BspHQ. +CI-PCR											
B	CN											
C	Sample 1 + CI-PCR											
D	...											
E												
F												
G												
H												

CP-BspHQ: positive control Bartonella spp/henselae/quintana; CI-PCR: control of PCR inhibition directly added in the Master Mix if not added before extraction; CN: Negative Control (water); Sample: Unknown DNA sample.

Protocol of qPCR

- 1) Homogenize the Master Mix and centrifuge briefly.
- 2) If the CI-PCR has not been added before extraction achieve the reaction mix as follows. Otherwise use directly 20µl of Master Mix.

The reaction mix should be carried out as follows, multiplying by the number of samples + standards (N) to test. Provide N + 3. Here the CI-PCR is added at the time of the PCR. It can also be added before extraction. It is available from Eurobio under reference EurobioPlex EBX - 002.

- 3) Distribute 20 µL of the reaction mix using a micropipette and sterile filter tips in each tube/well of microplate for real-time PCR.
- 4) Add 5 µL of DNA sample.
- 5) In parallel achieve :
 - a tube/well with 20 µL of the reaction mix + 5 µL of positive control CP-Barto (well A1)
 - a negative tube/well with 20 µL of master mix + 5 µL of water molecular biology (= H₂O: Well B1)
- 6) Close immediately with an adhesive film to avoid contamination.
- 7) Centrifuge briefly to collect the reaction mix at the bottom of the tubes or wells.
- 8) Set up the following program on qPCR instrument

Program	Temperature	Duration	Cycle(s)	
Activation of Taq polymerase	95°C	3 min	1	-
Amplification	95°C	10 sec	45	-
	60°C	30 sec		Fluorescence acquisition

Note 1 : Set heating and cooling rate by default, up to 20°C/sec, or 100%.

Note 2 : On LightCycler[®] 480 systems, two optical systems are available : only "System II" is compliant with use of the kit.

Note 3 : For the Applied Biosystems systems, select "NONE" in "PASSIVE REFERENCE".

Note 4 : On Rotorgene™, please calibrate the signal by clicking on "GAIN OPTIMIZATION".

VALIDATION OF THE EXPERIMENT

The threshold must be positioned right above the maximum of the negative control signal.

For the assay to be valid, the Ct values for the controls must be the following. Outside of these values, the experiment is not valid.

Table 3 :

Controls	CP-BsppHQ +CI-PCR	CN
Channel FAM	Ct ≤ 18	Ct not determined
Channel HEX	Ct ≤ 20	Ct not determined
Channel Texas Red	Ct ≤ 18	Ct not determined
Channel CY5	Ct ≤ 22	Ct not determined

CP-BsppHQ: positive control Bartonella spp/henselae/quintana; CI-PCR: control of PCR inhibition;
CN : negative control (water)

DATA ANALYSIS AND INTERPRETATION

For clinical samples, the following results are possible:

PCR Signal		Presence of bacteria	Test validity/comment
FAM or HEX or Texas Red	CY5	Bartonella spp Bartonella henselae Bartonella quintana	
+	+	Yes	valid
-	+	No	valid
+	-	Yes	Possible inhibition of PCR that doesn't prevent the detection of the bacteria; valid
-	-	?	Inhibition of PCR -dilute 5 x the sample; otherwise extract once again the sample

PERFORMANCES ANALYSIS

1. Studies of reproducibility:

Studies of reproducibility of the kit EurobioPlex Bartonella EBX-013 have been achieved on plasmid BsppHQ (positive control of Bartonella spp., Bartonella henselae and Bartonella quintana) on CFX96 (Biorad). The tables below indicate the average of the coefficient of variation (CV) depending on the concentration of the plasmid.

Reproducibility within an experiment (intra-experiment):

copies/ μ l	CV % CHANNEL FAM Bartonella spp	CV % CHANNEL HEX Bartonella henselae	CV % CHANNEL TEXAS RED Bartonella quintana
10 8	0.50	0,05	0,35
10 7	0.35	0,48	0,44
10 6	0.48	0,22	0,19
10 5	0.30	0,09	0,12
10 4	0.50	0,29	0,30
10 3	0.27	0,28	0,28
100	0.35	0,50	0,33
50	1.56	2,10	1,63
25	2.00	2,16	2,32
12.5	2.80	2,88	2,52
10			1,60
5	/	/	2,91
AVERAGE CV %	0,91	0,91	1,08


Reproducibility between experiments (inter-experiments):

copies/ μ l	CV % CHANNEL FAM Bartonella spp	CV % CHANNEL HEX Bartonella henselae	CV % CHANNEL TEXAS RED Bartonella quintana
10 8	8,91	3,96	2,80
10 7	6,36	1,99	4,04
10 6	3,77	1,32	4,62
10 5	3,11	1,30	3,49
10 4	2,58	1,09	3,27
10 3	0,19	0,48	4,86
10 2	1,38	2,08	3,89
10	1,24	2,40	4,00
5			4,26
AVERAGE CV %	3,44	1,83	3,91

2. Performance analysis :

Analytical sensitivity: 12,5 copies/ μ l for spp and henselae ; 5 copies/ μ l for quintana.
Linearity for quantification: from the detection limit to 10E+09 copies/ μ l.

Experiment conducted on CFX96 (BioRad) qPCR thermocycler:

	CHANNEL	Efficacy	Coefficient of correlation R 2	Slope
Barto spp	FAM	95.5 %	0.999	-3.436
Barto henselae	HEX	98.5 %	0.995	-3.358
Barto quintana	TEXAS RED	95.3 %	0.998	-3.440

3. Diagnostic evaluation: Clinical study

The tests were performed on 109 samples: 3 sera, 6 blood samples, 90 lymph node biopsies, 1 spleen biopsy, 1 extract from bone marrow, 7 cardiac biopsies, 1 sample of unknown origin). These samples were previously tested by the routine laboratory technique from the French National Reference Center (CNR) in Marseille and compared with the kit EurobioPlex Bartonella EBX-013 on CFX96 (Biorad).

The specificity and sensitivity of the EBX-013 kit have been evaluated by analysis of :

- 60 Bartonella spp POSITIVES
- 49 Bartonella spp NEGATIVES
- 45 Bartonella henselae POSITIVES/quintana NEGATIVES
- 64 Bartonella henselae NEGATIVES, 15 of them are the 15 quintana positives ; no cross-reactivity
- 15 Bartonella quintana POSITIVES/henselae NEGATIVES
- 94 Bartonella quintana NEGATIVES, 45 of them are the 45 henselae positives

		EBX-0.13 Quadriplex Bartonella EUROBIO								
		B spp			B henselae			B quintana		
		POS	NEG	TOTAL	POS	NEG	TOTAL	POS	NEG	TOTAL
Prestested CNR	POS	60	0	60	41	4	45	15	0	15
	NEG	2	47	49	0	64	64	0	94	94
	TOTAL	62	47	109	41	68	109	15	94	109

All samples were correctly analysed and all have a valid signal for internal control.

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WASTE DISPOSAL

Dispose of waste according to DASRI regulations.

SYMBOLS

	Reference
	Batch number
	Highest storage temperature
	Expiration Date
	Content sufficient for « N » reactions
	Keep protected from light
	Manufacturer
	CE labeled product
	In Vitro Diagnostic
	Instructions for use

eurobio *AbCys*

7, avenue de Scandinavie
ZA de Courtaboeuf
91953 Les Ulis Cedex
FRANCE