For Professional Use Only

**EuroClone**<sup>®</sup>  
diagnostica s.r.l.

## Duplica<sup>RealTime</sup> Mix & Match Factor V G1691A Kit

REV. EER038032\_IFU\_REV.00D\_ENITA

REF: EER038032- 32 tests

### Instructions For Use

#### INTENDED USE

**Duplica<sup>RealTime</sup> Mix & Match FACTOR V G1691A Kit** is an *in vitro* nucleic acid amplification test for the detection of the nucleotide substitution G1691A (allelic variant g.169549811C>T) of the blood coagulation Factor V gene in genomic DNA extracted from peripheral whole blood samples collected in EDTA.

EuroClone Diagnostica's **Duplica<sup>RealTime</sup> Mix & Match kits** belong to the CardioVascular Diseases panel – *CVD Panel* (EER037032, EER038032, EER039032 and EER040032) and share a common thermal profile.

#### INTRODUCTION

The resistance to Activated Protein C (APC), characterized by a poor response to anticoagulant activity of APC, a key enzyme in the downregulation of blood coagulation which cause a disposition for hypercoagulable state, increases the risk of thrombotic events. A number of clinical studies show a prevalence of APC resistance of 20-60% among patients with Venous ThromboEmbolic (VTE). The actual thrombotic risk is moderate with an odds ratio of 5-7 but its high prevalence makes it by far the most important risk factor know today, even higher than the sum of contributions from inherited deficiencies of Antithrombin, protein C and protein S. At least 90% of the cases with resistance to APC are explained by a nucleotide substitution in the gene for the coagulation Factor V (Guanine to Adenine, G1691A) causing the substitution of an Arginine to Glutamine at position 506 (FV Q506, often denoted as Factor V Leiden mutation), originating one of the three APC cleavage sites in activated Factor V. The mutation is inherited as an autosomic dominant trait and has a prevalence of 2% to 10% in the general caucasian population. Heterozygous carriers for Factor V Leiden have a 5 to 8 times increased risk of thrombosis and a 30 to 40 times for homozygous carriers. Recently it has been showed a 1.3-fold increased risk only for a coronary disease. Heterozygous carriers who take oral contraceptives are at a 15-fold increased risk of VTE, while carriers also heterozygous with Factor V Leiden have an approximate 20-fold higher risk and it seems to be higher when the coagulation Factor II G20210A transition is associated with the Factor V Leiden mutation.

#### PRINCIPLE OF THE TEST

**Duplica<sup>RealTime</sup> Mix & Match FACTOR V G1691A Kit** is designed to identify the G1691A substitution of the blood coagulation Factor V gene. Two reaction mixes are provided for the amplification:

- **AMPLIFICATION MIX**, containing Hot Start Taq DNA polymerase, nucleotides, MgCl<sub>2</sub> and buffer.
- **OLIGO MIX**, containing primers and fluorogenic probes.

**Duplica<sup>RealTime</sup> Mix & Match FACTOR V G1691A Kit** is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulation of PCR amplification products by fluorescent

DNA probes. In particular, the probe designed to detect the Wild Type allele carries the fluorophore FAM (6-carboxy-fluorescein) at the 5' end, while the probe detecting the Mutated allele, is labelled with the fluorophore HEX (hexachloro-fluorescein). Both the probes have a non fluorescent black quencher at the 3' end. If excited, the whole probe doesn't emit fluorescence, since the proximity of the quencher to the reporter prevents the emission of the fluorescence from the reporter (quenching effect).

#### REAGENTS PROVIDED

Each kit contains enough reagents to perform 32 tests when used in 4 analytical sessions with 5 **samples**, 1 **Wild Type control (Control 1, C1)**, 1 **Mutated control (Control 2, C2)** and one **Reaction Blank (BM)** each.

#### Kit Components

<i>Reagent</i>	<i>Color Code</i>	<i>Storage (range, °C)</i>	<i>Volume (µl)</i>	<i>Quantity (tubes)</i>
<b>Oligo Mix (OM)*</b>	Green Cap	-22÷-18	400	1
<b>Amplification Mix (AM)</b>	Blue Cap	-22÷-18	400	1
<b>Control 1 (C1, Wild Type)</b>	Red Cap	-22÷-18	>50	1
<b>Control 2 (C2, Mutated)</b>	Yellow Cap	-22÷-18	>50	1
<b>Reaction Blank (BM)</b>	White Cap	-22÷-18	>50	1

**\*protect the tube from direct light**

#### STORAGE AND HANDLING

All reagents must be stored at **-22÷-18°C** and can be used until the expiry date printed on the labels. Do not freeze and thaw the products more than six times.

#### MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to the specific handbook's section)
- Optical tubes or microplate for Real Time PCR
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 µl, 20-200 µl and 100-1000 µl)
- Disposable RNase/DNase-free tips with aerosol barriers
- Tube racks
- Desktop centrifuge
- PCR box
- Refrigerator
- Deep-freezer
- Thermalcycler for Real Time PCR

The kit has been optimized to be used on DX<sup>®</sup>/CFX (Bio-Rad), Rotor-Gene<sup>®</sup> Q (Qiagen) and Applied Biosystems<sup>®</sup> 7500 (Life Technologies<sup>™</sup>) thermalcyclers. Other makes and models should be fully tested and evaluated for optimal performance by the user before reporting results.

The equipments should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure an optimal performance.

#### PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood
- If required, Euroclone Diagnostica offers the necessary technical support for the correct use of the kit.
- Carefully read this Instruction For Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipments. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite

- Use powder-free gloves. Do not leave fingerprints on optical caps. Do not write on caps as this may cause an interference with fluorescent detection
- Materials containing or potentially-containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood
- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet (MSDS, available on request)
- The kit reagents, individual protective equipments, used materials, biological samples and test residuals must be disposed in accordance with local regulations
- Patient Drug treatment may interfere with the final result of the molecular biology analysis

## **OPERATING PROCEDURE**

### **a) DNA purification**

For genomic DNA purification EuroClone Diagnostica recommends to use:

- Duplicα Blood DNA kit (ref. EDI002250) and Duplicα NA Body Fluid kit (ref. EDI004200) with Duplicα® PREP Automatic Extractor (ref. EDI001) for automated purification.

For manual extraction/purification:

- Fasst DNA Releaser (ref. EMR057050) for **fresh** peripheral blood (i.e. stored for up to 24 hours at 2÷8°C).
- Spin DNA Purification Kit (ref. EMR061050) for **frozen** blood or stored at 2÷8°C for more than 24 hours.

Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

**Attention! Use blood samples in EDTA anti-coagulant solution.**

### **b) Thermalcycler Setup**

Refer to the specific handbook of the equipment used to set the thermal profile indicated in the **Thermal Profile Table**. We recommend to switch on the instrument and to set the thermal profile before preparing the reaction mix.

#### **DX®/CFX platform**

- Select *Create a new Experiment*
- In the *Protocol* section select *Create New* to set the thermal profile
- In *Plate Editor*: set FAM (Wild Type genotype, Allele 1) and HEX (Mutated genotype, Allele 2) as *Fluorophores* and define the sample name

#### **Applied Biosystems® 7500 platform**

- Select *Create a new Experiment*
- In the *Experiment Properties* select *Quantitation-Standard curve* as experiment type and *Taqman® Reagents* as reagent type
- Select the Run Mode (*7500* or *7500 Fast*) and the Ramp Speed (*Standard*)
- In the *Setup* menu select *Run Method* to set the thermal profile
- In the *Setup* menu select *Plate Setup-Define Targets and Samples* to select Samples and assign VIC target to MUT (Mutated genotype) and FAM target to WT (Wild-Type genotype). Leave the default value (*NFQ-MGB*) as quencher type
- In *Define Samples* insert the sample name
- Select the used wells in *Plate Setup-Assign Targets and Samples*
- Select *None* as passive reference

#### **Rotor-Gene® Q platform**

- Start the software and on the box *New Run* select *Advanced*
- Select a *new template* in *Empty Run* or a pre-existing one
- Select the Rotor Type of your instrument and then *Next*
- Type 25 µl in the reaction volume and then *Next*
- Select *Edit Profile* and set up the correct Thermal Profile as indicated in the Table below
- Select *Gain Optimisation* and then flag the option *Perform Optimisation before 1<sup>st</sup> acquisition*
- On *Channel Settings* select the green/yellow fluorophores and tube position "1" to perform the optimization. Then close the window and select *Next* and *Start Run*

### Thermal Profile Table (common for all the supported platforms)

Time	Temperature	Cycles
10 min	95°C	1
20 sec	95°C	35
30 sec	60°C	Fluorescence Acquisition

#### c) Preparation of PCR mix

The total reaction volume is **25 µl**.

For each experiment prepare a PCR mix for the **2 controls (C1 and C2)**, **1 Reaction Blank (BM)** and **n+1 samples**. The reagents of the PCR mix have to be mixed as indicated in the table below:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Extracted DNA	5

**The PCR mix has to be freshly prepared every time.**

After its preparation, aliquot **20 µl of Master Mix** in the tubes or in the microplates well for PCR then add in each tube/well **5 µl** (100-250 ng/reaction) from the **extracted DNA** or **control DNA**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermalcycler.

#### d) ANALYSIS and INTERPRETATION of the RESULTS

##### Real Time PCR curves analysis

Refer to the instrument specific user guide to visualize the amplification plots for the entire plate/rotor. Detailed analysis of raw data depends on the Real Time PCR instrument used. Baseline noise levels should either be set automatically or at predefined cycles.

The fluorescence in each channel indicates the hybridisation of the allelic specific probes: **Channel 1** for **FAM/green= Wild Type allele probe** and **Channel 2** for **HEX/VIC/Yellow= Mutated allele probe**. If a sample shows a fluorescence in **FAM/Green**, the sample has the **Wild Type allele**. If a sample shows a fluorescence in **HEX/VIC/Yellow**, the sample has the **Mutated allele**.

Therefore, if only a **FAM/Green signal** is detected the sample is **Homozygous Wild Type**, whereas if only a **HEX/VIC/Yellow signal** is detected the sample is **Homozygous Mutated**. Finally, if both **FAM/Green** and **HEX/VIC/Yellow** are detected the sample is **Heterozygous**.

Condition in which no signal is detected indicates PCR inhibition. In this case, refer to troubleshooting.

Control 1 (Wild Type), Control 2 (Mutant) and reaction blank (BM) are provided in order to properly set the threshold line before samples analysis. After the run, the threshold line has to be set so that: **Control 1 is detected in FAM/Green and not detected in HEX/VIC/Yellow, whereas Control 2 is detected in HEX/VIC/Yellow and not detected in FAM/Green. The Reaction Blank must not be detected in any channel.**

**If these three conditions have been met, the run is valid and it's possible to analyse the data; otherwise the run is not valid. It's responsibility of the user to validate the run.**

#### Results Interpretation Table

FAM/Green	HEX/VIC/Yellow	Results
Detected	Not detected	Wild Type
Not detected	Detected	Mutated
Detected	Detected	Heterozygous
Not detected	Not detected	Inhibition

#### Scatter Plot analysis


Additionally, the **Scatter Plot Analysis** can be used to visualize the results for **Rotor-Gene®Q** and **DX®/CFX**

platforms. In this analysis each sample is identified by a spot representing the end-point fluorescence values for the two dyes.

### Rotor-Gene® Q platform


- In the *Analysis* menu select *Other*
- Select *Scatter Graph Analysis*, highlight both *Green* and *Yellow* channels and press *Show* to display the scatter plot. The scatter plot is divided in four areas corresponding to samples classified as Wild Type, Mutant, Heterozygous.

In the *Scatter Analysis graph* window that appears select with the mouse the area corresponding to a specific Genotype and name it as show in the figure below (*Green* → *Wild Type*; *Yellow* → *Mutated*; and *Green/Yellow* → *Heterozygous*):

### DX®/CFX platform

- Select *Ct Determination Mode* and *Regression* in the *Settings* menu
- In the *View/Edit Plate* menu set the reaction blank as *NTC*
- Select *Quantitation*. In the *Ct* table in the lower right part of the screen no amplification must be detected (*Ct* value=N/A for both fluorophores) for the **blank/NTC** sample. **Control 1 (C1)** must be detected (*Ct* value≠N/A) in FAM and non detected (*Ct* value=N/A) in HEX, whereas **Control 2 (C2)** detected in HEX and not detected in FAM. If these conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run by checking that these conditions have been met
- **Unknown (Unk)** samples must be detected in at least one fluorescence channel. Condition in which no signal is detected (*Ct* value=N/A for both fluorophores) indicates PCR inhibition and **samples presenting this condition must be excluded from the following Allelic Discrimination analysis**. For troubleshooting refer to the point 1 and 2 of the Troubleshooting section.

- In the *Allelic Discrimination* menu select *RFU* and flag the *Normalize Data* in the *Display Mode* to display the relative Scatter Plot
  - For both Fluorophores set manually the Threshold to 0,15
- The Scatter Plot displays three main quadrants in which each spot represents the two fluorescent signals for a single sample/control:

The corresponding genotype of the samples is reported in a table displayed close to the Scatter Plot:

<b>Bottom Right Quadrant</b>	Wild Type Samples (WT- <i>Allele 1</i> )
<b>Upper Left Quadrant</b>	Mutated Samples (MUT- <i>Allele 2</i> )
<b>Upper Right Quadrant</b>	Heterozygote Samples (HET- <i>both Alleles</i> )

## **TROUBLESHOOTING**

### **Problem 1: Weak or no signal of unknown samples.**

1. The PCR was inhibited:
  - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.
2. Pipetting error due to omitted reagents or samples:
  - Repeat the analysis starting from the PCR.
3. Deterioration of dyes and/or primers. The reagents storage conditions didn't comply with the instructions:
  - Check storage conditions.
4. Very low starting amount and/or low purity of genomic DNA. Improper DNA extraction:
  - Repeat the DNA purification.
5. Wrong channel/filter was chosen. The PCR conditions didn't comply with the instructions:
  - Check the PCR conditions and select the fluorescence channels reported in the protocol for the Unknown Sample detection.

### **Problem 2: Weak or no signal of the Wild Type Control C1 and/or Mutated Control C2.**

1. The PCR conditions didn't comply with the instructions:
  - Check the amplification protocol and select the fluorescence channel reported in the manual.
2. Deterioration of dyes and/or primers. The reagents storage conditions didn't comply with the instructions:
  - Check storage conditions.

### **Problem 3: Any signal with Reaction Blank BM.**

1. Contamination during PCR preparation procedure. All samples results are INVALID:
  - Decontaminate the working area and all instruments.
  - Pipette the controls C1 and C2 at last.
  - Repeat the PCR preparation using a new set of reagents.

### **Problem 4: HEX signal in Wild Type Control C1 and/or FAM signal in Mutated Control C2.**

1. Contamination during PCR preparation procedure. All samples results are INVALID:
  - Decontaminate the working area and all instruments.
  - Use tips with aerosol barriers only. Change tips between tubes/wells.
  - Repeat the PCR preparation using a new set of reagents.

### **Problem 5: Fluorescence intensity varies.**

1. The PCR Master Mix is not well prepared:
  - Carefully repeat the PCR preparation procedure.
2. Air bubbles trapped in the PCR tubes:
  - Check the presence of air bubbles before starting a new run.

### **Problem 6: Absence of any fluorescent signal.**

1. Verify the performance of the thermalcycler:
  - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions didn't comply with the instructions:
  - Check storage conditions.
  - Check the expiry date of the kit.

### **Problem 7: The thermalcycler gives an error message.**

1. Refer to the Real Time PCR instrument user manual or contact the local technical support of the Real Time PCR instrument company.

### **Problem 8: The kit reagents left out of the storage range temperature.**

1. These reagents must be stored as indicated for a proper execution of the test. The performance of the product is not guaranteed if the reagents have not been properly stored.

**"THE PURCHASE OF THIS PRODUCT GRANTS THE PURCHASER RIGHTS UNDER CERTAIN ROCHE PATENTS TO USE IT SOLELY FOR PROVIDING HUMAN IN VITRO DIAGNOSTIC SERVICES. NO GENERAL PATENT OR OTHER LICENSE OF ANY KIND OTHER THAN THIS SPECIFIC RIGHT OF USE FROM PURCHASE IS GRANTED HEREBY"**

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EuroClone Diagnostica Srl