
# REF M41 MICROGEN™ C. DIFFICILE

## INTENDED USE

Microgen™ C. difficile is a rapid latex agglutination test intended for confirmatory identification of *Clostridium difficile* cultured on selective solid media from faecal samples from patients with suspected pseudomembranous colitis, antibiotic-associated diarrhoea and post-operative diarrhoea. The kit is intended for professional laboratory use only.

## PRINCIPLE OF THE TEST

Latex particles are coated with rabbit IgG antibodies specific for *C. difficile* cell wall antigens. When the sensitised latex particles are mixed with a suspension of *C. difficile* colonies, a sensitive and specific immunochemical reaction takes place causing the finely dispersed latex particles to agglutinate rapidly into aggregates that are easily visible to the unaided eye.

CONT	KIT PRESENTATION
------	------------------

REAG	TEST	M41a Test Latex Reagent:	2.5mL
------	------	--------------------------	-------

Latex particles coated with rabbit antibodies to *C. difficile* antigens. Preserved with 0.099% sodium azide. (Blue cap)

CONTROL	+	M41b Positive Control:	0.5mL
---------	---	------------------------	-------

Suspension of inactivated *C. difficile* antigens reactive with Test Latex Reagent. Preserved with 0.099% sodium azide. (Black cap)

NaCl	0.85%	M40 0.85% Isotonic Saline:	5mL
------	-------	----------------------------	-----

Preserved with 0.099% sodium azide. (White cap)

Instructions for Use  
Disposable agglutination slides  
Disposable mixing sticks

## Additional Requirements:

- Bacteriological loops
- *C. difficile* selective agar plates (CCFA or CCEY)

## WARNINGS AND PRECAUTIONS

### Safety:

1. The reagents supplied in this kit are for *in vitro* diagnostic use only
2. Sodium azide, which is used as a preservative in the kit reagents can react with lead or copper plumbing to form potentially explosive metal azides. Dispose by flushing with a large volume of water to prevent azide build-up.
3. Appropriate precautions should be taken when handling or disposing of potential pathogens. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30 minutes. Liquid waste containing acid must be neutralised before treatment.

4. The positive control has been inactivated during the manufacturing process. However, it should be handled as though potentially infectious.

## Procedural:

1. Microgen™ C. difficile should be used according to the kit instructions.
2. Allow all reagents to reach room temperature before use.
3. Do not dilute any of the kit reagents
4. Do not intermix reagents from different batches of kits.
5. Do not freeze any of the kit reagents
6. Do not allow the latex reagent dropper to touch the positive control or bacterial samples.
7. Ensure the agglutination slide is clean and dry prior to use.
8. The kit should not be used if the latex reagent fails to agglutinate with the positive control, or if the latex reagent agglutinates in isotonic saline only. Replace with a new kit.

## STORAGE AND SHELF LIFE

Microgen™ C. difficile should be stored at 2-8°C when not in use. The kit should not be used after the expiry date printed on the carton label.

## SPECIMENS

Faecal samples should be cultured in enrichment broth at 37°C for 18-24 hours, then sub-cultured on to selective agar plates (CCFA or CCEY). Plates should be incubated anaerobically at 37°C for 24-48 hours.

Colonies with morphology resembling *C. difficile* are removed for testing with Microgen™ C. difficile (see below)

## PROCEDURE

### Quality Control:

The following checks should be performed each time the kit is used to confirm that the reagents are functioning correctly:

1. **Reagent Control**  
Add 1 drop of Microgen™ C. difficile Latex Reagent to 1 drop of isotonic saline in the same circle on an agglutination slide. Mix with a mixing stick and observe for agglutination. No agglutination should be seen. If this control shows agglutination, at least one of the reagents is contaminated and they should be discarded.
2. **Positive Control**  
Gently mix the Positive Control by inverting several times. Place 1 drop on a circle of an agglutination slide. Add 1 drop of Microgen™ C. difficile Latex Reagent to the same circle and mix. Agglutination should be visible within 2 minutes. If no agglutination is seen the reagents should be discarded.

### Test Procedure:

1. Dispense 1 drop of isotonic saline on to 1 circle of a clean, dry Microgen™ agglutination slide.
2. Using an inoculating loop, remove a suspected *C. difficile* colony from the selective agar plate. Only select colonies whose morphology resembles that of *C. difficile*. Emulsify the colony in the drop of saline on the test card to produce a heavy, smooth suspension.

- Observe the suspension for any agglutination or clumping which would indicate auto-agglutination. If the suspension remains smooth, proceed to the next step. **If auto-agglutination is seen, the organism cannot be tested using Microgen™ C. difficile. Alternative test methods should be used.**
- Gently mix the Test Latex Reagent by inverting the vial several times. Add 1 drop to the colony suspension on the slide. **Do not allow the dropper to touch the organism suspension.**
- Mix the latex reagent and organism suspension together with a clean mixing stick for 30 seconds. Continue mixing by rocking the slide.
- Examine for agglutination after 2 minutes from initial mixing of latex and sample.
- After reading, discard the used slides and mixing sticks into suitable disinfectant.

## INTERPRETATION

Agglutination within 2 minutes is a positive result and indicates the presence of *C. difficile*.

No agglutination within 2 minutes is a negative result.

## LIMITATIONS OF USE

- Results should be interpreted by the clinician in the context of all available clinical and laboratory information. The isolation of *C. difficile* does not constitute a diagnosis of pseudomembranous colitis or antibiotic-associated diarrhoea.
- Identification of *C. difficile* using Microgen™ should be performed on selective cultures as this increases the isolation rate.
- Other *Clostridium spp* may occasionally grow on CCFA selective agar. However, they can be distinguished from *C. difficile* by culture on CCEY agar. Some will not grow on this medium whilst others do not display the characteristics of *C. difficile*.
- Culture-derived suspensions which auto-agglutinate cannot be tested by Microgen™ *C. difficile*. Alternative methods should be used.

## PERFORMANCE CHARACTERISTICS

Microgen™ *C. difficile* has been evaluated as a culture confirmation test at both an independent UK microbiology laboratory and in-house. In total, 137 bacterial isolates were cultured on selective agar plates and colonies tested by Microgen™ *C. difficile* and a well-established commercially available test.

		Microgen™ <i>C. difficile</i>		
		+ve	-ve	Total
Commercial test	+ve	85*	0	85
	-ve	0	52**	52
Total		85	52	137

Sensitivity: 85/85 = 100%  
 Specificity: 52/52 = 100%  
 Diagnostic Efficiency: 137/137 = 100%

\*Of the 85 isolates in this group, 18 were cross-reactants in both tests. However, 16 of these either will not grow on *C. difficile* selective medium or the colonies do not resemble *C. difficile*. The remaining two isolates (both *C. glycolicum*) will grow slightly but do not exhibit colony fluorescence which is a characteristic of the majority of *C. difficile* strains.

\*\* 2 of these isolates were classified as *C. difficile* (serogroups A9, A10). One of these (serogroup A10) exhibited slightly irregular colony morphology. The remaining 50 isolates comprised a wide variety of bacterial species including 5 *Clostridium spp*. Most of these isolates either do not grow on *C. difficile* selective agar or exhibit atypical colony morphology.

Overall, the results obtained with Microgen™ *C. difficile* correlate closely with those obtained using the established commercial product. Although a number of organisms have the potential to cause false

positive reactions in both tests, they either do not grow in *C. difficile* selective culture media or their colony morphologies are not typical of *C. difficile*.

## REPRODUCIBILITY

**Intra-batch reproducibility** was established by testing one batch of product on three separate occasions using a different operator for each occasion. Sensitivity was tested using serial dilutions of reference and kit control antigens and specificity was confirmed using a QC organism panel. No substantial differences were seen between the results obtained by the three operators.

**Inter-batch reproducibility** was examined by testing the sensitivity and specificity of three batches of product using reference and kit control antigens and the QC organism panel. No differences in sensitivity or specificity were seen between the three batches of product.

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Microgen Bioproducts Ltd