
# Midi Parasep® SF Faecal Parasite Concentrator

<b>CS</b>	Midi Parasep® SF Koncentrátor parazitů ve stolici
<b>DE</b>	Midi Parasep® SF Konzentratör für Stuhlparasiten
<b>ES</b>	Midi Parasep® SF Concentrador de parásitos fecales
<b>FR</b>	Midi Parasep® SF Concentrateurs de Parasites Fécaux
<b>HR</b>	Midi Parasep® SF Koncentratori crijevnih parazita bez otapala
<b>IT</b>	Midi Parasep® SF Concentratore di parassiti fecali
<b>NL</b>	Midi Parasep® SF Fecale Parasieten Concentrator
<b>PL</b>	Midi Parasep® SF System do zagęszczania kału przy analizie parazytów
<b>PT</b>	Midi Parasep® SF Concentrador de parasitas fecais
<b>SI</b>	Midi Parasep® SF Koncentrator parazitov v blatu

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See label for storage conditions and expiry date. Please adhere to the following guidelines when handling Midi Parasep® SF. To avoid cross contamination the Midi Parasep® SF device should remain closed at all times except when introducing the sample or when retrieving the final concentrated sample for examination.

### Sample Preparation

If using prefilled Midi Parasep® SF, start at 1D

- 1A Unscrew lid.
- 1B Add 8.0ml of fixative.
- 1C Add one drop of surfactant (eg: Apacor Triton X solution) to the chamber.
- 1D Introduce a pea sized faecal sample to the fixative. Mix in thoroughly with the Midi Parasep® SF spoon. If the sample is hard, break it up with the end of the spoon.

### Emulsification

- 2 Seal Midi Parasep® SF by screwing in the filter/sedimentation cone unit. Vortex or shake to emulsify with the sedimentation cone pointing upwards.

### Centrifugation

- 3 Invert Midi Parasep® SF and centrifuge at 400g for 2 minutes.  
(J. Clin. Microbiol. doi:10.1128/JCM.00838-15)  
Midi Parasep® SF fits all 50ml centrifuge buckets.

NOTE: To calculate the required RPM for any centrifuge

$$\text{RPM} = \sqrt{\frac{g}{1.12r}} \times 1000$$

RPM Rotor Speed in revs/min

g centrifugal force (max. 1000g)

r radius, horizontal distance between sedimentation cone tip and spindle centre measured in mm

### Examination

- 4A Unscrew and discard the filter and mixing tube.
- 4B Pour off all the liquid above the sediment.
- 4C Pipette one drop of sediment onto a slide and cover with coverslip.  
Alternatively, follow laboratory SOP for slide preparation.

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**ALCORFIX™ SAFETY DATA SHEET**

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**SECTION 1 IDENTIFICATION OF THE SUBSTANCE/ MIXTURE AND THE COMPANY/UNDERTAKING****1.1 Product Identifier: AlcorFix™****108810, 108885, 108886, 108887, 108889, 149995, 248200, 249200, 249300, 249420****1.2 Relevant identified uses of the substance or mixture and uses advised against:** Solution for fixation/conservation of biological samples.**1.3 Details of the supplier of the Safety Data Sheet:**

Apacor Limited, Unit 5 Sapphire Centre, Fishponds Road,  
Wokingham, Berkshire, RG41 2QL, England  
+44 (0) 118 979 5566  
[technical@apacor.com](mailto:technical@apacor.com)

**1.4 Emergency telephone number:**

+44 (0)118 979 5566  
(Monday-Friday 0900-1700 excluding UK Public Holidays)

**SECTION 2 HAZARDS IDENTIFICATION****2.1 Classification of the substance or mixture**

2.1.1 Classification according to Regulation (EC) No 1272/2008 [CLP]:

Acute toxicity, Oral (Category 4), H302

Acute toxicity, Inhalation (gas) (Category 4), H332

Serious eye damage (Category 1), H318

Hazardous to the aquatic environment (Category 2), H411

Flammable liquids (Category 2), H225

See Section 16 for the full text of H-Statements mentioned in this Section.

**2.2 Label elements**

Labelling according to Regulation (EC) No 1272/2008 [CLP]

Pictogram

Signal word

Danger

**Hazard statement(s)**

H225 – Highly flammable liquid and vapour

H302 - Harmful if swallowed

H318 - Causes serious eye damage

H332 - Harmful if inhaled

H411 - Toxic to aquatic life with long lasting effects

**Precautionary statements:**

P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P280 - Wear protective gloves/protective clothing/eye protection/face protection.

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 - Immediately call a POISON CENTER or doctor/physician.

P370 + P378 - In case of fire: Use dry sand, carbon dioxide (CO<sub>2</sub>), water spray, dry chemical or alcohol resistant foam to extinguish.

**2.3 Other hazards**

None.

**SECTION 3 COMPOSITION/INFORMATION ON INGREDIENTS****3.1 Substances****3.2 Mixtures****Hazardous ingredients according to Regulation (EC) No 1272/2008****Component: Ethanol**

CAS No: 64-17-5

EC No: 200-578-6

Index No: 603-002-00-5

Classification: Flam. Liq. 2 (H225)

Concentration: 25%

**Component: Zinc sulphate**

CAS No: 7733-02-0

EC No: 231-793-3

Index No: 030-006-00-9

Classification: Acute Tox. 4 (H302), Eye Dam. 1 (H318), Aquatic Acute 1 (H400), Aquatic Chronic 1 (H410)

Concentration: 7.9%

**Component: Acetic Acid**

CAS No: 64-19-7

EC No: 200-580-7

Index No: 607-002-00-6

Classification: Skin Corr. 1A (H314), Flam. Liq. 3 (H226)

Concentration: 4.8%

**Component: Isopropanol**

CAS No: 67-63-0

EC No: 200-661-7

Index No: 603-117-00-0

Classification: Eye Irrit. 2 (H319), STOT SE 3 (H336), Flam. Liq. 2 (H225)

Concentration: 1%

**Component: Methyl Alcohol**

CAS No: 67-56-1

EC No: 200-659-6

Index No: 603-001-00-X

Classification: Acute Tox. 3 (H301), Acute Tox. 3 (H311), Acute Tox. 3. (H331), STOT SE 1 (H370), Flam. Liq. 2 (H225)

Concentration: 1%

**3.3 Other Information**

Additional non-hazardous ingredients:

Polyvinyl alcohol (minimum 1g/l)

DI water

**SECTION 4 FIRST AID MEASURES****4.1 Description of first aid measures**

Consult a physician. Show this safety data sheet to the doctor in attendance.

**In case of eye contact:** Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

## ALCORFIX™ SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**In case of skin contact:** Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes.

**If swallowed:** Clean mouth with water and drink afterwards plenty of water.

**If inhaled:** Move to fresh air.

#### 4.2 Most important symptoms and effects, both acute and delayed

No information available.

#### 4.3 Indication of any immediate medical attention and special treatment needed

**Notes to physician:** Treat symptomatically.

### SECTION 5 FIRE FIGHTING MEASURES

#### 5.1 Extinguishing media

Suitable extinguishing media (use media appropriate to the circumstances and environment): dry sand, carbon dioxide (CO<sub>2</sub>), water spray, alcohol-resistant foam, dry chemical.

#### 5.2 Special hazards arising from the substance or mixture

No information available

#### 5.3 Advice for firefighters

As in any fire, wear self-contained breathing apparatus, MSHA/NIOSH (approved or equivalent) and full protective gear.

### SECTION 6 ACCIDENTAL RELEASE MEASURES

#### 6.1 Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation, especially in confined areas.

#### 6.2 Environmental precautions

Should not be released into the environment. Prevent product from entering drains. Prevent further leakage or spillage if safe to do so.

#### 6.3 Methods and material for containment and cleaning up

Absorb spill with inert material (eg dry sand or earth), then place in a chemical waste container. After cleaning, flush away traces with water.

#### 6.4 Reference to other sections

For disposal, see Section 13.

### SECTION 7 HANDLING AND STORAGE

#### 7.1 Precautions for safe handling

Do not breathe vapours or spray mist. Ensure that ventilation is adequate before using this product. Avoid contact with skin and eyes. Take necessary personal protective precautions before using this product. Keep away from heat and flame. Take precautionary measures against static discharges.

#### 7.2 Conditions for safe storage, including any incompatibilities

**Technical measures/precautions:** Store in a tightly closed container. Store in a cool, dry, well-ventilated area away from incompatible substances.

**Incompatible products:** Avoid strong bases. Oxidizing agent.

#### 7.3 Specific end use(s)

No other specific end uses(s) are specified apart from those listed in Section 1.2.

### SECTION 8 EXPOSURE CONTROLS/PERSONAL PROTECTION

#### 8.1 Control parameters

Component	Ethanol 64-17-5	Zinc sulphate 7733-02-0	Acetic Acid 64-19-7	Isopropanol 67-63-0	Methyl Alcohol 67-56-1
UK	STEL: 3000 ppm STEL: 5760 mg/m <sup>3</sup> TWA: 1000 ppm TWA: 1920 mg/m <sup>3</sup>			STEL: 500 ppm STEL: 1250 mg/m <sup>3</sup> TWA: 400 ppm TWA: 999 mg/m <sup>3</sup>	STEL: 250 ppm STEL: 333 mg/m <sup>3</sup> TWA: 200 ppm TWA: 266 mg/m <sup>3</sup> Skin
France	TWA: 1000 ppm TWA: 1900 mg/m <sup>3</sup> STEL: 5000 ppm STEL: 9500 mg/m <sup>3</sup>		STEL: 10 ppm STEL: 25 mg/m <sup>3</sup>	STEL: 400 ppm STEL: 980 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 260 mg/m <sup>3</sup> STEL: 1000 ppm STEL: 1300 mg/m <sup>3</sup>
Spain	STEL: 1000 ppm STEL: 1910 mg/m <sup>3</sup>		STEL: 15 ppm STEL: 37 mg/m <sup>3</sup> TWA: 10 ppm TWA: 25 mg/m <sup>3</sup>	STEL: 400 ppm STEL: 1000 mg/m <sup>3</sup> TWA: 200 ppm TWA: 500 mg/m <sup>3</sup>	S* TWA: 200 ppm TWA: 266 mg/m <sup>3</sup>
Germany	TWA: 500 ppm TWA: 960 mg/m <sup>3</sup> Ceiling / Peak: 1000 ppm Ceiling / Peak: 1920 mg/m <sup>3</sup> Skin	TWA: 0.1 mg/m <sup>3</sup> TWA: 2 mg/m <sup>3</sup> Ceiling / Peak: 0.4 mg/m <sup>3</sup> Ceiling / Peak: 4 mg/m <sup>3</sup>	TWA: 10 ppm TWA: 25 mg/m <sup>3</sup> Ceiling / Peak: 20 ppm Ceiling / Peak: 50 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 500 mg/m <sup>3</sup> Ceiling / Peak: 400 ppm Ceiling / Peak: 1000 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 270 mg/m <sup>3</sup> Ceiling / Peak: 800 ppm Ceiling / Peak: 1080 mg/m <sup>3</sup> Skin
Italy					TWA: 200 ppm TWA: 260 mg/m <sup>3</sup> Skin
Portugal	TWA: 1000 ppm		STEL: 15 ppm TWA: 10 ppm TWA: 25 mg/m <sup>3</sup>	STEL: 400 ppm TWA: 200 ppm	STEL: 250 ppm TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
The Netherlands	Skin STEL: 1900 mg/m <sup>3</sup> TWA: 260 mg/m <sup>3</sup>				Skin TWA: 133 mg/m <sup>3</sup> TWA: 100 ppm
Finland	TWA: 1000 ppm TWA: 1900 mg/m <sup>3</sup> STEL: 1300 ppm STEL: 2500 mg/m <sup>3</sup>		TWA: 5 ppm TWA: 13 mg/m <sup>3</sup> STEL: 10 ppm STEL: 25 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 500 mg/m <sup>3</sup> STEL: 250 ppm STEL: 620 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 270 mg/m <sup>3</sup> STEL: 250 ppm STEL: 330 mg/m <sup>3</sup> Skin
Denmark	TWA: 1000 ppm TWA: 1900 mg/m <sup>3</sup>		TWA: 10 ppm TWA: 25 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 490 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 260 mg/m <sup>3</sup> Skin
Austria	STEL 2000 ppm STEL 3800 mg/m <sup>3</sup> TWA: 1000 ppm TWA: 1900 mg/m <sup>3</sup>		STEL 20 ppm STEL 50 mg/m <sup>3</sup> TWA: 10 ppm TWA: 25 mg/m <sup>3</sup>	STEL 800 ppm STEL 2000 mg/m <sup>3</sup> TWA: 200 ppm TWA: 500 mg/m <sup>3</sup>	Skin STEL 800 ppm STEL 1040 mg/m <sup>3</sup> TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
Switzerland	STEL: 1000 ppm STEL: 1920 mg/m <sup>3</sup> TWA: 500 ppm TWA: 960 mg/m <sup>3</sup>	STEL: 4 mg/m <sup>3</sup> TWA: 0.1 mg/m <sup>3</sup> TWA: 2 mg/m <sup>3</sup>	STEL: 20 ppm STEL: 50 mg/m <sup>3</sup> TWA: 10 ppm TWA: 25 mg/m <sup>3</sup>	STEL: 400 ppm STEL: 1000 mg/m <sup>3</sup> TWA: 200 ppm TWA: 500 mg/m <sup>3</sup>	Skin STEL: 800 ppm STEL: 1040 mg/m <sup>3</sup> TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
Poland	TWA: 1900 mg/m <sup>3</sup>		STEL: 30 mg/m <sup>3</sup> TWA: 15 mg/m <sup>3</sup>	STEL: 1200 mg/m <sup>3</sup> TWA: 900 mg/m <sup>3</sup>	TWA: 300 mg/m <sup>3</sup> TWA: 100 mg/m <sup>3</sup>
Norway	TWA: 500 ppm TWA: 950 mg/m <sup>3</sup> STEL: 500 ppm STEL: 950 mg/m <sup>3</sup>		TWA: 10 ppm TWA: 25 mg/m <sup>3</sup> STEL: 20 ppm STEL: 37.5 mg/m <sup>3</sup>	TWA: 100 ppm TWA: 245 mg/m <sup>3</sup> STEL: 150 ppm STEL: 306.25 mg/m <sup>3</sup>	TWA: 100 ppm TWA: 130 mg/m <sup>3</sup> Skin STEL: 150 ppm STEL: 162.5 mg/m <sup>3</sup>
Ireland	STEL: 1000 ppm		TWA: 10 ppm TWA: 25 mg/m <sup>3</sup> STEL: 15 ppm STEL: 37 mg/m <sup>3</sup>	TWA: 200 ppm STEL: 400 ppm Skin	TWA: 200 ppm TWA: 260 mg/m <sup>3</sup> STEL: 600 ppm STEL: 780 mg/m <sup>3</sup> Skin
European Union				TWA 10 ppm TWA 25 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 260 mg/m <sup>3</sup> Skin

Derived No Effect Level (DNEL) No information available

Predicted No Effect Concentration (PNEC) No information available

#### 8.2 Exposure controls

**Engineering measures:** Ensure adequate ventilation, especially in confined areas.

#### Personal protective equipment

**Respiratory protection:** No special protective equipment required.

**Hand protection:** Wear appropriate protective gloves.

**Eye protection:** Wear tightly fitting safety goggles or safety glasses with side-shields.

**Skin and body protection:** Protective clothing to protect exposed skin.

**Hygiene measures:** Handle in accordance with good industrial hygiene and safety practice.

**Environmental exposure controls:** No information available.

## ALCORFIX™ SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES****9.1 Information on basic physical and chemical properties**

- a) Appearance:** clear liquid
- b) Odour:** pungent
- c) Odour threshold:** no information available
- d) pH:** no information available
- e) Melting point / freezing point:** no information available
- f) Initial boiling point / boiling range:** 84°C
- g) Flash point:** 16°C
- h) Evaporation rate:** no information available
- i) Flammability (solid, gas):** no information available
- j) Upper/lower flammability or explosive limits:** no information available
- k) Vapour pressure:** no information available
- l) Vapour density:** no information available
- m) Relative density:** no information available
- n) Solubility (ies):** soluble in water
- o) Partition coefficient: n-octanol/water:** no information available
- p) Auto-ignition temperature:** no information available
- q) Decomposition temperature:** no information available
- r) Viscosity:** no information available
- s) Explosive properties:** no information available
- t) Oxidising properties:** no information available
- 9.2 Other information:** no information available

**SECTION 10 STABILITY AND REACTIVITY****10.1 Reactivity****10.2 Chemical stability**

Stable under normal conditions.

**10.3 Possibility of hazardous reactions**

No information available.

**10.4 Conditions to avoid**

Heat, flames and sparks.

**10.5 Incompatible materials**

No particular materials.

**10.6 Hazardous decomposition products**

Under normal use – none.

**SECTION 11 TOXICOLOGICAL INFORMATION****11.1 Information on toxicological effects****Acute toxicity:**

Product: based on known/supplied information, does not present an acute toxicity hazard.

Inhalation: no data available.

Eye contact: no data available.

Skin contact: no data available.

Ingestion: no data available.

&lt; 60.3% of the mixture consists of ingredients of unknown toxicity.

The following values are calculated based on GHS document chapter 3.1.

Oral	1,363.00mg/kg
Dermal	5,158.00mg/kg
Inhalation:	
Gas	4,263.00mg/l
Mist	20.90mg/l
Vapour	829.22mg/l

Chemical Name	LD50 Oral	LD50 Dermal	LC50 Inhalation
Ethanol	7060mg/kg (Rat)		124.7mg/L (Rat) 4 h
Zinc sulphate	500mg/kg (Rat)		
Acetic acid	3310mg/kg (Rat)	1060mg/kg (Rabbit)	11.4mg/L (Rat) 4 h
Methyl alcohol	6200mg/kg (Rat)	15800mg/kg (Rabbit)	22500 ppm (Rat) 8 h
Isopropanol	1870mg/kg (Rat)	4059mg/kg (Rabbit)	72600mg/m <sup>3</sup> (Rat) 4 h

**Skin corrosion/irritation:** no data available**Serious eye damage/eye irritation:** no data available**Respiratory or skin sensitisation:** no data available**Germ cell mutagenicity:** no data available**Carcinogenicity:** no data available**Reproductive toxicity:** no data available**Specific target organ toxicity - single exposure:** no data available**Specific target organ toxicity - repeated exposure:** no data available**Aspiration hazard:** no data available**SECTION 12 ECOLOGICAL INFORMATION****12.1 Toxicity**

Toxic to aquatic life with long lasting effects.

Chemical Name	Toxicity to Algae	Toxicity to Fish	Toxicity to Daphnia and other aquatic invertebrates
Ethanol		12.0 - 16.0: 96 h Oncorhynchus mykiss m/L LC50 static 100: 96 h Pimephales promelas mg/L LC50 static 13400 - 15100: 96 h Pimephales promelas mg/L LC50 flow-through	9268 - 14221: 48 h Daphnia magna mg/L LC50 2: 48 h Daphnia magna mg/L EC50 Static 10800: 24 h Daphnia magna mg/L EC50
Zinc sulphate	0.056: 72 h Pseudokirchneriella subcapitata mg/L EC50 static 64.8: 72 h Chlorella vulgaris mg/L EC50 2.4: 96 h Chlorella vulgaris mg/L EC50	0.162: 96 h Oncorhynchus mykiss mg/L LC50 flow-through 0.03 - 0.05: 96 h Oncorhynchus mykiss mg/L LC50 semi-static 0.34 - 0.93: 96 h Oncorhynchus mykiss mg/L LC50 static 0.218 - 0.42: 96 h Pimephales promelas mg/L LC50 flow-through 0.06: 96 h Pimephales promelas mg/L LC50 static 0.23 - 0.48: 96 h Pimephales promelas mg/L LC50 0.168 - 0.25: 96 h Pimephales promelas mg/L LC50 semi-static 0.15: 96 h Cyprinus carpio mg/L LC50 semi-static 16.85 - 27.18: 96 h Cyprinus carpio mg/L LC50 static 3 - 4.6: 96 h Lepomis macrochirus mg/L LC50 flow-through 3.55 - 6.32: 96 h Lepomis macrochirus mg/L LC50 static 0.63: 96 h Poecilia reticulata mg/L LC50 mg/L LC50 49.23 - 64.16: 96 h Poecilia reticulata mg/L LC50 semi-static 0.48 - 1.72: 96 h Poecilia reticulata mg/L LC50 static	0.75: 48 h Daphnia magna mg/L EC50 0.538 - 0.908: 48 h Daphnia magna mg/L EC50 Static
Acetic acid		79: 96 h Pimephales promelas mg/L LC50 static 75: 96 h Lepomis macrochirus mg/L LC50 static	65: 48 h Daphnia magna mg/L EC50 Static 47: 24 h Daphnia magna mg/L EC50
Isopropanol	1000: 96 h Desmodesmus subspicatus mg/L EC50 1000: 72 h Desmodesmus subspicatus mg/L EC50	9640: 96 h Pimephales promelas mg/L LC50 flow-through 11130: 96 h Pimephales promelas mg/L LC50 static 1400000: 96 h Lepomis macrochirus µg/L LC50	13299: 48 h Daphnia magna mg/L EC50
Methyl alcohol		28200: 96 h Pimephales promelas mg/L LC50 flow-through 100: 96 h Pimephales promelas mg/L LC50 static 19500 - 20700: 96 h Oncorhynchus mykiss mg/L LC50 flow-through 18 - 20: 96 h Oncorhynchus mykiss mL/L LC50 static 13500 - 17600: 96 h Lepomis macrochirus mg/L LC50 flow-through	

**ALCORFIX™ SAFETY DATA SHEET**

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**12.2 Persistence and degradability**

No information available.

**12.3 Bioaccumulative potential**

No information available.

Chemical Name	log Pow
Ethanol	-0.32
Acetic acid	-0.31
Isopropanol	0.05
Methyl alcohol	-0.77

**12.4 Mobility in soil**

No information available.

**12.5 Results of PBT and vPvB assessment**

No information available

**12.6 Other adverse effects**

No information available.

**12.7 Additional information**

No information available.

**SECTION 13 DISPOSAL CONSIDERATIONS****13.1 Waste treatment methods**

**Waste from residues / unused products:** In accordance with local and national regulations. Should not be released into the environment.

**Contaminated packaging:** Empty containers should be disposed of at an approved waste handling site for recycling or disposal.

**SECTION 14 TRANSPORT INFORMATION****14.1 UN number:** UN2924

**14.2 UN proper shipping name:** Flammable Liquid, Corrosive, n.o.s. (Ethanol, Acetic Acid)

**14.3 Transport hazard class(es):** 3, Subsidiary Class: 8

**14.4 Packing group:** II

**14.5 Environmental hazards****14.6 Special precautions for user****14.7 Transport in bulk according to Annex II of****MARPOL73/78 and the IBC Code**

Not intended to be transported in bulk.

Note: Per 49 CFR – when shipping 30ml or less per inner packaging and the gross weight does not exceed 64lbs, use the 173.4 small quantity exception.

**SECTION 15 REGULATORY INFORMATION****15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

Chemical Name	French RG number
Ethanol	RG 84
Isopropanol	RG 84
Methyl alcohol	RG 84

TSCA	Complies
EINECS/ELINCS	-
DSL/NDSL	-
PICCS	-
ENCS	-
IECSC	-
AICS	-
KECL	-

**Legend**

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory  
EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances  
DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List  
PICCS - Philippines Inventory of Chemicals and Chemical Substances  
ENCS - Japan Existing and New Chemical Substances  
IECSC - China Inventory of Existing Chemical Substances  
AICS - Australian Inventory of Chemical Substances  
KECL - Korean Existing and Evaluated Chemical Substances

**15.2 Chemical Safety Assessment**

No information available.

**SECTION 16 OTHER INFORMATION****Full text of H-Statements referred to under sections 2 and 3**

H225 - Highly flammable liquid and vapour  
H226 - Flammable liquid and vapour  
H301 - Toxic if swallowed  
H302 - Harmful if swallowed  
H311 - Toxic in contact with skin  
H314 - Causes severe skin burns and eye damage  
H318 - Causes serious eye damage  
H319 - Causes serious eye irritation  
H331 - Toxic if inhaled  
H332 - Harmful if inhaled  
H336 - May cause drowsiness or dizziness  
H370 - Causes damage to organs.  
H400 - Very toxic to aquatic life  
H410 - Very toxic to aquatic life with long lasting effects  
H411 - Toxic to aquatic life with long lasting effects

Amended sections are indicated by a line in the border.

The information supplied in this SDS is correct to the best of our knowledge. We do not accept any liability for loss, injury or damage, which may result from its use.

## 10% FORMALIN SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

### **SECTION 1 IDENTIFICATION OF THE SUBSTANCE/ MIXTURE AND THE COMPANY/UNDERTAKING**

#### **1.1 Product Identifier: 10% Formalin**

145200, 145300, 145400, 145420, 145700, 145800, 145900, 1460, 146200, 146300, 146400, 108900, 108910, 148980, 148998, 149910, 151000

#### **1.2 Relevant identified uses of the substance or mixture and uses advised against:** laboratory chemical (in vitro diagnostic)

#### **1.3 Details of the supplier of the Safety Data Sheet:**

Apacor Limited, Unit 5 Sapphire Centre, Fishponds Road, Wokingham, Berkshire, RG41 2QL, England

+44 (0) 118 979 5566

[technical@apacor.com](mailto:technical@apacor.com)

#### **1.4 Emergency telephone number:**

+44 (0)118 979 5566

(Monday-Friday 0900-1700 excluding UK Public Holidays)

### **SECTION 2 HAZARDS IDENTIFICATION**

#### **2.1 Classification of the substance or mixture**

Classification according to Regulation (EC) No 1272/2008 [CLP]:

Acute toxicity, Oral (Category 4), H302

Skin sensitisation (Category 1), H317

Acute toxicity, Inhalation (Category 4), H332

Germ cell mutagenicity (Category 2), H341

Carcinogenicity (Category 1B), H350

See Section 16 for the full text of H-Statements mentioned in this Section.

#### **2.2 Label elements**

Labelling according to Regulation (EC) No 1272/2008 [CLP]

Danger

#### **Pictogram**

#### **Signal word**

#### **Hazard statement(s)**

H302 Harmful if swallowed

H317 May cause an allergic skin reaction.

H332 Harmful if inhaled.

H341 Suspected of causing genetic defects

H350 May cause cancer

Contains Formaldehyde.

#### **Precautionary statements:**

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301 + P312 + P330 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Rinse mouth.

P308 + P313 IF exposed or concerned: Get medical advice/attention.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

See Section 16 for the full text of H-Statements mentioned in this Section.

#### **2.3 Other hazards**

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

### **SECTION 3 COMPOSITION/INFORMATION ON INGREDIENTS**

#### **3.2 Mixtures**

#### **Hazardous ingredients according to Regulation (EC) No 1272/2008**

Component: **Formaldehyde**

CAS No: 50-00-0

EC No: 200-001-8

Index No: 605-001-00-5

Classification: Acute Tox. 3 (H301 + H311 + H331), Skin Corr. 1B (H314), Skin Sens. 1 (H317), Muta. 2 (H341), Carc. 1B (H350)

Concentration: < 5%

Component: **Methanol**

CAS No: 67-56-1

EC No: 200-659-6

Index No: 603-001-00-x

Registration No: 01-2119433307-44-xxxx

Classification: Flam. Liq. 2 (H225); Acute Tox. 3 (H301 + H311 + H331); STOT SE 1 H370

Concentration: < 1%

See Section 16 for the full text of H-Statements mentioned in this Section.

### **SECTION 4 FIRST AID MEASURES**

#### **4.1 Description of first aid measures**

Consult a physician. Show this safety data sheet to the doctor in attendance.

**If inhaled:** If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

**In case of skin contact:** Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

**In case of eye contact:** Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

**If swallowed:** Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

#### **4.2 Most important symptoms and effects, both acute and delayed**

The most important known symptoms and effects are described in the labelling (Section 2.2) and/or Section 11.

#### **4.3 Indication of any immediate medical attention and special treatment needed**

No data available.

### **SECTION 5 FIRE FIGHTING MEASURES**

#### **5.1 Extinguishing media**

Suitable extinguishing media: Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

# 10% FORMALIN SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

## 5.2 Special hazards arising from the substance or mixture

Carbon oxides

## 5.3 Advice for firefighters

Wear self-contained breathing apparatus and full protective gear.

## SECTION 6 ACCIDENTAL RELEASE MEASURES

### 6.1 Personal precautions, protective equipment and emergency procedures

Use personal protective equipment. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Remove all sources of ignition. Evacuate personnel to safe areas. Beware of vapours accumulating to form explosive concentrations. Vapours can accumulate in low areas. For personal protection see Section 8.

### 6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

### 6.3 Methods and material for containment and cleaning up

Contain spillage, and then collect and place in container for disposal according to local regulations (see Section 13). Keep in suitable, closed containers for disposal.

### 6.4 Reference to other sections

For disposal, see Section 13.

## SECTION 7 HANDLING AND STORAGE

### 7.1 Precautions for safe handling

Avoid contact with skin and eyes. Avoid inhalation of vapour or mist. Keep away from sources of ignition—no smoking. Take measures to prevent the build-up of electrostatic charge. For precautions see Section 2.2.

### 7.2 Conditions for safe storage, including any incompatibilities

Store in a cool place. Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

### 7.3 Specific end use(s)

No other specific uses are specified apart from those listed in Section 1.2.

## SECTION 8 EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1 Control parameters

Exposure limits: this product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

	<b>Formaldehyde 50-00-0</b>	<b>Methanol 67-56-1</b>
<b>Austria</b>	STEL: 0.5 ppm STEL: 0.6 mg/m <sup>3</sup> TWA: 0.5 ppm TWA: 0.6 mg/m <sup>3</sup>	STEL: 800 ppm STEL: 1040 mg/m <sup>3</sup> TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
<b>Belgium</b>	STEL: 0.3 ppm STEL: 0.38 mg/m <sup>3</sup>	STEL: 250 ppm STEL: 333 mg/m <sup>3</sup> TWA: 200 ppm TWA: 266 mg/m <sup>3</sup>
<b>Denmark</b>	STEL: 0.3 ppm STEL: 0.4 mg/m <sup>3</sup> TWA: 0.3 ppm TWA: 0.4 mg/m <sup>3</sup>	STEL: 400 ppm STEL: 520 mg/m <sup>3</sup> TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>

	<b>Formaldehyde 50-00-0</b>	<b>Methanol 67-56-1</b>
<b>France</b>	TWA: 0.5 ppm STEL: 1 ppm	STEL: 1000 ppm STEL: 1300 mg/m <sup>3</sup> TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
<b>Germany</b>	STEL: 0.6 ppm STEL: 0.74 mg/m <sup>3</sup> TWA: 0.3 ppm TWA: 0.37 mg/m <sup>3</sup>	STEL: 800 ppm STEL: 1080 mg/m <sup>3</sup> TWA: 200 ppm TWA: 270 mg/m <sup>3</sup>
<b>Ireland</b>	STEL: 2 ppm STEL: 2.5 mg/m <sup>3</sup> TWA: 2 ppm TWA: 2.5 mg/m <sup>3</sup>	TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
<b>Italy</b>		TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
<b>Poland</b>	STEL: 1 mg/m <sup>3</sup> TWA: 0.5 mg/m <sup>3</sup>	STEL: 300 mg/m <sup>3</sup> TWA: 100 mg/m <sup>3</sup>
<b>Portugal</b>	STEL: 0.3 ppm	STEL: 250 ppm TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>
<b>Spain</b>	STEL: 0.3 ppm STEL: 0.37 mg/m <sup>3</sup>	STEL: 250 ppm STEL: 333 mg/m <sup>3</sup> TWA: 200 ppm TWA: 266 mg/m <sup>3</sup>
<b>Sweden</b>	STEL: 0.6 ppm STEL: 0.74 mg/m <sup>3</sup> TWA: 0.3 ppm TWA: 0.37 mg/m <sup>3</sup>	STEL: 250 ppm STEL: 350 mg/m <sup>3</sup> TWA: 200 ppm TWA: 250 mg/m <sup>3</sup>
<b>The Netherlands</b>	STEL: 0.5 mg/m <sup>3</sup> TWA: 0.15 mg/m <sup>3</sup>	TWA: 133 mg/m <sup>3</sup>
<b>UK</b>	STEL: 2 ppm STEL: 2.5 mg/m <sup>3</sup> TWA: 2 ppm TWA: 2.5 mg/m <sup>3</sup>	STEL: 250 ppm STEL: 333 mg/m <sup>3</sup> TWA: 200 ppm TWA: 266 mg/m <sup>3</sup>

### 8.2 Exposure controls

#### 8.2.1 Appropriate engineering controls

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

#### 8.2.2 Personal protective equipment

(a) Eye/face protection: Tightly fitting safety goggles. Faceshield (8-inch minimum). Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

(b) Skin protection: Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves should satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

(c) Body Protection: Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

(d) Respiratory protection: Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

## 10% FORMALIN SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

### 8.2.3 Environmental exposure controls

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

## SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

a) **Appearance** Form: liquid

b) **Odour** no data available

c) **Odour threshold** no data available

d) **pH** no data available

e) **Melting point / freezing point** no data available

f) **Initial boiling point and boiling range** 100°C at 1.013 hPa

g) **Flash point** 85°C

h) **Evaporation rate** no data available

i) **Flammability (solid, gas)** no data available

j) **Upper/lower flammability or explosive limits**

Upper 70% (V), Lower 7% (V)

k) **Vapour pressure** 53hPa at 39°C

l) **Vapour density** no data available

m) **Relative density** 1.080g/cm<sup>3</sup>

n) **Solubility (ies)** completely miscible

o) **Partition coefficient: n-octanol/water** no data available

p) **Auto-ignition temperature** no data available

q) **Decomposition temperature** no data available

r) **Viscosity** no data available

s) **Explosive properties** no data available

t) **Oxidising properties** no data available

### 9.2 Other information

No data available.

## SECTION 10 STABILITY AND REACTIVITY

### 10.1 Reactivity

No data available.

### 10.2 Chemical stability

Stable under recommended storage conditions.

### 10.3 Possibility of hazardous reactions

No data available.

### 10.4 Conditions to avoid

Heat, flames and sparks.

### 10.5 Incompatible materials

No materials to be mentioned in particular.

### 10.6 Hazardous decomposition products

Carbon oxides.

## SECTION 11 TOXICOLOGICAL INFORMATION

### 11.1 Information on toxicological effects

**Acute toxicity:** no data available

**Skin corrosion/irritation:** no data available

**Serious eye damage/eye irritation:** no data available

**Respiratory or skin sensitisation:** no data available

**Germ cell mutagenicity:** no data available

**Carcinogenicity:** IARC: 1 - Group 1: Carcinogenic to humans (Formaldehyde)

**Reproductive toxicity:** no data available

**Specific target organ toxicity - single exposure:** no data available

**Specific target organ toxicity - repeated exposure:** no data available

**Aspiration hazard:** no data available

### Additional Information

#### Chemical Name

Formaldehyde	LD50 oral 600mg/kg (Rat) LD50 dermal 270mg/kg (Rabbit) LC50 inhalation 0.578mg/L (Rat) 4 h
Methanol	LD50 oral - rat - 5628mg/kg LC50 inhalation - rat - 4h – 83.2mg/l/4h

## SECTION 12 ECOLOGICAL INFORMATION

### 12.1 Toxicity

Ecotoxicity effects: contains no substances known to be hazardous to the environment or not degradable in waste water treatment plants.

#### Toxicity to Fish

Formaldehyde	0.032 - 0.226: 96 h Oncorhynchus mykiss mL/L LC50 flow-through 100- 136: 96 h Oncorhynchus mykiss mg/L LC50 static 1510: 96 h Lepomis macrochirus µg/L LC50 static 22.6 - 25.7: 96 h Pimephales promelas mg/L LC50 flow-through 23.2 - 29.7: 96 h Pimephales promelas mg/L LC50 static 41: 96 h Brachydanio rerio mg/L LC50 static
Methanol	LC50 - Pimephales promelas – 28200mg / L 96h

#### Toxicity to Daphnia and other Aquatic Invertebrates

Formaldehyde	11.3 - 18: 48 h Daphnia magna mg/L EC50 Static 2: 48 h Daphnia magna mg/L LC50
Methanol	EC50 - Daphnia magna - >10000mg/l

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

#### Chemical Name log Pow

Formaldehyde	0.35
Methanol	-0.77

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

### 12.6 Other adverse effects

No data available.

### 12.7 Additional information

None.

**10% FORMALIN SAFETY DATA SHEET**

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**SECTION 13 DISPOSAL CONSIDERATIONS****13.1 Waste treatment methods**

**Product:** Dispose of in accordance with all federal, state, and local regulations. This combustible material may be burned in a chemical incinerator equipped with an afterburner and scrubber. Offer surplus and non-recyclable solutions to a licensed disposal company.

**Contaminated packaging:** Dispose of as unused product.

Amended sections are indicated by a line in the border.

The information supplied in this SDS is correct to the best of our knowledge. We do not accept any liability for loss, injury or damage, which may result from its use.

**SECTION 14 TRANSPORT INFORMATION**

**IATA/DOT/IMDG/TDG:** Not regulated.

**14.1 UN number:** -

**14.2 UN proper shipping name:** -

**14.3 Transport hazard class(es):** -

**14.4 Packing group:** -

**14.5 Environmental hazards:** -

**14.6 Special precautions for user:** -

**SECTION 15 REGULATORY INFORMATION****15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

No data available.

**15.2 Chemical Safety Assessment**

No chemical safety assessment has been carried out for this product.

**SECTION 16 OTHER INFORMATION****Full text of H-Statements referred to in Sections 2 and 3**

H225 Highly flammable liquid and vapour.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

H317 May cause an allergic skin reaction.

H332 Harmful if inhaled.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

H370 Causes damage to organs.

Acute Tox. Acute toxicity

Carc. Carcinogenicity

Flam. Liq. Flammable liquids

Muta. Germ cell mutagenicity.

Skin Corr. Skin corrosion

Skin Sens. Skin sensitisation

STOT SE Specific target organ toxicity - single exposure

**SAF SAFETY DATA SHEET**

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**SECTION 1 IDENTIFICATION OF THE SUBSTANCE/ MIXTURE AND THE COMPANY/UNDERTAKING****1.1 Product Identifier: SAF**

(Sodium Acetate-Acetic Acid-Formalin Solution)

145500, 1461, 1461A, 146500, 146501, 108920, 148960,  
148965, 149920, 149960, 153000, 249400**1.2 Relevant identified uses of the substance or mixture and uses advised against:** laboratory chemical (in vitro diagnostic)**1.3 Details of the supplier of the Safety Data Sheet:**Apacor Limited, Unit 5 Sapphire Centre, Fishponds Road,  
Wokingham, Berkshire, RG41 2QL, England

+44 (0) 118 979 5566

[technical@apacor.com](mailto:technical@apacor.com)**1.4 Emergency telephone number:**

+44 (0)118 979 5566

(Monday-Friday 0900-1700 excluding UK Public Holidays)

**SECTION 2 HAZARDS IDENTIFICATION****2.1 Classification of the substance or mixture**

Classification according to Regulation (EC) No 1272/2008

[CLP]:

Acute toxicity, Oral (Category 4), H302

Skin sensitisation (Category 1), H317

Acute toxicity, Inhalation (Category 4), H332

Germ Cell Mutagenicity (Category 2), H341

Carcinogenicity (Category 1B), H350

See Section 16 for the full text of H-Statements mentioned in this Section.

**2.2 Label elements**

Labelling according to Regulation (EC) No 1272/2008 [CLP]

**Pictogram****Signal word****Danger****Hazard statement(s)**

H302 Harmful if swallowed.

H317 May cause an allergic skin reaction.

H332 Harmful if inhaled.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

Contains Formaldehyde

**Precautionary statements:**

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician

P308 + P313 - IF exposed or concerned: Get medical advice/ attention

**2.3 Other hazards**

No data available.

**SECTION 3 COMPOSITION/INFORMATION ON INGREDIENTS****3.2 Mixtures****Hazardous ingredients according to Regulation (EC) No 1272/2008**Component: **Formaldehyde**

CAS No: 50-00-0

EC No: 200-001-8

Index No: 605-001-00-5

Classification: Acute Tox. 3 (H301 + H311 + H331); Skin Corr. 1B (H314); Skin Sens. 1 (H317); Muta. 2 (H341); Carc. 1B (H350); Concentration: &lt; 5%

Component: **Methanol**

CAS No: 67-56-1

EC No: 200-659-6

Index No: 603-001-00-x

Registration No: 01-2119433307-44-xxxx

Classification: Flam. Liq. 2 (H225); Acute Tox 3 (H301 + H311 + H331); STOT SE 1 (H370)

Concentration: &lt; 1%

Component: **Acetic Acid**

CAS No: 64-19-7

EC No: 200-580-7

Index No: -

Registration No: -

Classification: Skin Corr. 1A (H314) ; Flam. Liq 3 (H226)

Concentration: ≤ 2%

**SECTION 4 FIRST AID MEASURES****4.1 Description of first aid measures**

Consult a physician. Show this safety data sheet to the doctor in attendance.

**If inhaled:** If breathed in, move person into fresh air. If not breathing, give artificial respiration.**In case of skin contact:** Wash off immediately with soap and plenty of water for at least 15 minutes while removing all contaminated clothes and shoes.**In case of eye contact:** Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.**If swallowed:** Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.**4.2 Most important symptoms and effects, both acute and delayed**

The most important known symptoms and effects are described in the labelling (Section 2.2) and/or in Section 11.

**4.3 Indication of any immediate medical attention and special treatment needed**

Notes to physician: treat symptomatically.

**SECTION 5 FIRE FIGHTING MEASURES****5.1 Extinguishing media**

Suitable extinguishing media: Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

**5.2 Special hazards arising from the substance or mixture**

Carbon oxides.

**SAF SAFETY DATA SHEET**

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**5.3 Advice for firefighters**

Wear self-contained breathing apparatus and full protective gear.

**SECTION 6 ACCIDENTAL RELEASE MEASURES****6.1 Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Remove all sources of ignition. Evacuate personnel to safe areas. Beware of vapours accumulating to form explosive concentrations.

Vapours can accumulate in low areas. For personal protection see Section 8.

**6.2 Environmental precautions**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

**6.3 Methods and material for containment and cleaning up**

Contain spillage and place in container for disposal according to local regulations (see Section 13). Keep in suitable, closed containers for disposal.

**6.4 Reference to other sections**

For disposal, see Section 13.

**SECTION 7 HANDLING AND STORAGE****7.1 Precautions for safe handling**

Avoid contact with skin and eyes. Avoid inhalation of vapour or mist. Keep away from sources of ignition—no smoking. Take measures to prevent the build-up of electrostatic charge. For precautions, see Section 2.2.

**7.2 Conditions for safe storage, including any incompatibilities**

Store in a cool place. Keep container tightly closed in a dry and well-ventilated place. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

**7.3 Specific end use(s)**

No other specific uses are specified apart from those listed in Section 1.2.

**SECTION 8 EXPOSURE CONTROLS/PERSONAL PROTECTION****8.1 Control parameters**

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

	<b>Formaldehyde 50-00-0</b>	<b>Methanol 67-56-1</b>	<b>Acetic Acid 64-19-7</b>
<b>Austria</b>	STEL: 0.5 ppm	STEL: 800 ppm	STEL: 20 ppm
	STEL: 0.6 mg/m <sup>3</sup>	STEL: 1040 mg/m <sup>3</sup>	STEL: 50 mg/m <sup>3</sup>
	TWA: 0.5 ppm	TWA: 200 ppm	TWA: 10 ppm
	TWA: 0.6 mg/m <sup>3</sup>	TWA: 260 mg/m <sup>3</sup>	TWA: 25 mg/m <sup>3</sup>
<b>Belgium</b>	STEL: 0.3 ppm	STEL: 250 ppm	STEL: 15 ppm
	STEL: 0.38 mg/m <sup>3</sup>	STEL: 333 mg/m <sup>3</sup>	STEL: 38 mg/m <sup>3</sup>
		TWA: 200 ppm	TWA: 10 ppm
		TWA: 266 mg/m <sup>3</sup>	TWA: 25 mg/m <sup>3</sup>
<b>Denmark</b>	STEL: 0.3 ppm	STEL: 400 ppm	STEL: 20 ppm
	STEL: 0.4 mg/m <sup>3</sup>	STEL: 520 mg/m <sup>3</sup>	STEL: 50 mg/m <sup>3</sup>
	TWA: 0.3 ppm	TWA: 200 ppm	TWA: 10 ppm
	TWA: 0.4 mg/m <sup>3</sup>	TWA: 260 mg/m <sup>3</sup>	TWA: 25 mg/m <sup>3</sup>
<b>France</b>	TWA: 0.5 ppm	STEL: 1000 ppm	STEL: 10 ppm
	STEL: 1 ppm	STEL: 1300 mg/m <sup>3</sup>	STEL: 25 mg/m <sup>3</sup>
		TWA: 200 ppm	
		TWA: 260 mg/m <sup>3</sup>	

	<b>Formaldehyde 50-00-0</b>	<b>Methanol 67-56-1</b>	<b>Acetic Acid 64-19-7</b>
<b>Germany</b>	STEL: 0.6 ppm	STEL: 800 ppm	STEL: 20 ppm
	STEL: 0.74 mg/m <sup>3</sup>	STEL: 1080 mg/m <sup>3</sup>	STEL: 50 mg/m <sup>3</sup>
	TWA: 0.3 ppm	TWA: 200 ppm	TWA: 10 ppm
	TWA: 0.37 mg/m <sup>3</sup>	TWA: 270 mg/m <sup>3</sup>	TWA: 25 mg/m <sup>3</sup>
<b>Ireland</b>	STEL: 2 ppm	TWA: 200 ppm	STEL: 15 ppm
	STEL: 2.5 mg/m <sup>3</sup>	TWA: 260 mg/m <sup>3</sup>	STEL: 37 mg/m <sup>3</sup>
	TWA: 2 ppm		TWA: 10 ppm
	TWA: 2.5 mg/m <sup>3</sup>		TWA: 25 mg/m <sup>3</sup>
<b>Italy</b>		TWA: 200 ppm	TWA: 10 ppm
		TWA: 260 mg/m <sup>3</sup>	TWA: 25 mg/m <sup>3</sup>
<b>Poland</b>	STEL: 1 mg/m <sup>3</sup>	STEL: 300 mg/m <sup>3</sup>	STEL: 30 mg/m <sup>3</sup>
	TWA: 0.5 mg/m <sup>3</sup>	TWA: 100 mg/m <sup>3</sup>	TWA: 15 mg/m <sup>3</sup>
<b>Portugal</b>	STEL: 0.3 ppm	STEL: 250 ppm	STEL: 15 ppm
		TWA: 200 ppm	TWA: 10 ppm
		TWA: 260 mg/m <sup>3</sup>	TWA: 25 mg/m <sup>3</sup>
<b>Spain</b>	STEL: 0.3 ppm	STEL: 250 ppm	STEL: 15 ppm
	STEL: 0.37 mg/m <sup>3</sup>	STEL: 333 mg/m <sup>3</sup>	STEL: 37 mg/m <sup>3</sup>
		TWA: 200 ppm	TWA: 10 ppm
		TWA: 266 mg/m <sup>3</sup>	TWA: 25 mg/m <sup>3</sup>
<b>Sweden</b>	STEL: 0.6 ppm	STEL: 250 ppm	STEL: 10 ppm
	STEL: 0.74 mg/m <sup>3</sup>	STEL: 350 mg/m <sup>3</sup>	STEL: 25 mg/m <sup>3</sup>
	TWA: 0.3 ppm	TWA: 200 ppm	TWA: 5 ppm
	TWA: 0.37 mg/m <sup>3</sup>	TWA: 250 mg/m <sup>3</sup>	TWA: 13 mg/m <sup>3</sup>
<b>The Netherlands</b>	STEL: 0.5 mg/m <sup>3</sup>	TWA: 133 mg/m <sup>3</sup>	
	TWA: 0.15 mg/m <sup>3</sup>		
<b>UK</b>	STEL: 2 ppm	STEL: 250 ppm	
	STEL: 2.5 mg/m <sup>3</sup>	STEL: 333 mg/m <sup>3</sup>	
	TWA: 2 ppm	TWA: 200 ppm	
	TWA: 2.5 mg/m <sup>3</sup>	TWA: 266 mg/m <sup>3</sup>	

**8.2 Exposure controls****8.2.1 Appropriate engineering controls**

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

**8.2.2 Personal protective equipment**

(a) Eye/face protection: Tightly fitting safety goggles. Faceshield (8-inch minimum). Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

(b) Skin protection: Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves should satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

(c) Body Protection: Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

(d) Respiratory protection: Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

**8.2.3 Environmental exposure controls**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

**SAF SAFETY DATA SHEET**

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES**

- 9.1 Information on basic physical and chemical properties**
- a) **Appearance** aqueous solution Form: colourless liquid
  - b) **Odour** characteristic
  - c) **Odour threshold** no data available
  - d) **pH** no data available
  - e) **Melting point / freezing point** no data available
  - f) **Initial boiling point and boiling range** 102°C
  - g) **Flash point** >105°C
  - h) **Evaporation rate** no data available
  - i) **Flammability (solid, gas)** no data available
  - j) **Upper/lower flammability or explosive limits** no data available
  - k) **Vapour pressure** no data available
  - l) **Vapour density** >1
  - m) **Relative density** 1.071
  - n) **Solubility (ies)** Soluble in water
  - o) **Partition coefficient: n-octanol/water** no data available
  - p) **Auto-ignition temperature** no data available
  - q) **Decomposition temperature** no data available
  - r) **Viscosity** no data available
  - s) **Explosive properties** no data available
  - t) **Oxidising properties** no data available

**9.2 Other information**

No data available.

**SECTION 10 STABILITY AND REACTIVITY****10.1 Reactivity**

No data available.

**10.2 Chemical stability**

Stable under recommended storage conditions.

**10.3 Possibility of hazardous reactions**

No data available.

**10.4 Conditions to avoid**

Heat, flames and sparks.

**10.5 Incompatible materials**

No materials to be mentioned in particular.

**10.6 Hazardous decomposition products**

Carbon oxides.

**SECTION 11 TOXICOLOGICAL INFORMATION****11.1 Information on toxicological effects****Acute toxicity:** no data available**Skin corrosion/irritation:** no data available**Serious eye damage/eye irritation:** no data available**Respiratory or skin sensitisation:** no data available**Germ cell mutagenicity:** no data available**Carcinogenicity:** IARC: 1 - Group 1: Carcinogenic to humans  
(Formaldehyde)**Reproductive toxicity:** no data available**Specific target organ toxicity - single exposure:** no data available**Specific target organ toxicity - repeated exposure:** no data available**Aspiration hazard:** no data available**Additional Information****Chemical Name**

Formaldehyde	LD50 oral 600 mg/kg (Rat) LD50 dermal 270 mg/kg (Rabbit) LC50 inhalation 0.578 mg/L (Rat) 4 h
Methanol	LD50 oral - rat - 5628 mg / kg LC50 inhalation - rat - 4h – 83.2 mg/l/4h
Acetic Acid	LD50 oral 3310 mg/kg (Rat) LD50 dermal 1060 mg/kg (Rabbit) LC50 inhalation 11.4 mg/L (Rat) 4 h

**SECTION 12 ECOLOGICAL INFORMATION****12.1 Toxicity**

Ecotoxicity effects: contains no substances known to be hazardous to the environment or not degradable in waste water treatment plants.

**Toxicity to Fish**

Formaldehyde	0.032 - 0.226: 96 h Oncorhynchus mykiss mL/L LC50 flow-through 100- 136: 96 h Oncorhynchus mykiss mg/L LC50 static 1510: 96 h Lepomis macrochirus µg/L LC50 static 22.6 - 25.7: 96 h Pimephales promelas mg/L LC50 flow-through 23.2 - 29.7: 96 h Pimephales promelas mg/L LC50 static 41: 96 h Brachydanio rerio mg/L LC50 static
Methanol	LC50 - Pimephales promelas – 28200 mg / L 96h
Acetic Acid	75: 96 h Lepomis macrochirus mg/L LC50 static 79: 96 h Pimephales promelas mg/L LC50 static

**Toxicity to Daphnia and other Aquatic Invertebrates**

Formaldehyde	11.3 - 18: 48 h Daphnia magna mg/L EC50 Static 2: 48 h Daphnia magna mg/L LC50
Methanol	EC50 - Daphnia magna - >10000 mg/l
Acetic Acid	47: 24 h Daphnia magna mg/L EC50 65: 48 h Daphnia magna mg/L EC50 Static

**12.2 Persistence and degradability**

No data available.

**12.3 Bioaccumulative potential**

No data available.

**Chemical Name**

Chemical Name	log Pow
Formaldehyde	0.35
Methanol	-0.77
Acetic Acid	0

**12.4 Mobility in soil**

No data available.

**12.5 Results of PBT and vPvB assessment**

No data available.

**12.6 Other adverse effects**

No data available.

**12.7 Additional information**

None.

**SAF SAFETY DATA SHEET**

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

**SECTION 13 DISPOSAL CONSIDERATIONS****13.1 Waste treatment methods**

**Product:** Dispose of waste in accordance with all federal, state, and local regulations.

**Contaminated packaging:** Dispose of as unused product.

Amended sections are indicated by a line in the border.

The information supplied in this SDS is correct to the best of our knowledge. We do not accept any liability for loss, injury or damage, which may result from its use.

**SECTION 14 TRANSPORT INFORMATION**

**IATA/DOT/IMDG/TDG:** Not regulated.

**14.1 UN number:** -

**14.2 UN proper shipping name:** -

**14.3 Transport hazard class(es):** -

**14.4 Packing group:** -

**14.5 Environmental hazards:** -

**14.6 Special precautions for user:** -

**SECTION 15 REGULATORY INFORMATION****15.1 Safety, health and environmental regulations/  
legislation specific for the substance or mixture**

No data available.

**15.2 Chemical Safety Assessment**

No chemical safety assessment has been carried out for this product.

**SECTION 16 OTHER INFORMATION****Full text of H-Statements referred to in Sections 2 and 3**

H225 Highly flammable liquid and vapour.

H226 Flammable liquid and vapour.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

H317 May cause an allergic skin reaction.

H332 Harmful if inhaled.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

H370 Causes damage to organs.

Acute Tox. Acute toxicity

Carc. Carcinogenicity

Flam. Liq. Flammable liquids

Muta. Germ Cell Mutagenicity

Skin Corr. Skin corrosion

Skin Sens. Skin sensitisation

STOT SE Specific target organ toxicity - single exposure

## TRITON X SOLUTION SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

### **SECTION 1 IDENTIFICATION OF THE SUBSTANCE/ MIXTURE AND THE COMPANY/UNDERTAKING**

#### **1.1 Product Identifier: 1472, 172018**

#### **TRITON X Solution**

Used at concentration of ≤0.1% in: 145300, 145400, 145420, 145500, 145800, 145900, 146300, 146400, 146500, 146501, 108900, 108910, 108920, 108935, 148960, 148965, 148980, 148998, 149910, 149920, 149960, 151000, 153000, 249400

#### **1.2 Relevant identified uses of the substance or mixture and uses advised against:** for laboratory use (in vitro diagnostic).

#### **1.3 Details of the supplier of the Safety Data Sheet:**

Apacor Limited, Unit 5 Sapphire Centre, Fishponds Road, Wokingham, Berkshire, RG41 2QL, England  
+44 (0) 118 979 5566  
[technical@apacor.com](mailto:technical@apacor.com)

#### **1.4 Emergency telephone number:**

+44 (0)118 979 5566  
(Monday-Friday 0900-1700 excluding UK Public Holidays)

### **SECTION 2 HAZARDS IDENTIFICATION**

#### **2.1 Classification of the substance or mixture**

Classification according to Regulation (EC) No 1272/2008 [CLP]:

Serious eye damage (Category1), H318

See Section 16 for the full text of H-Statements mentioned in this Section.

#### **2.2 Label elements**

Labelling according to Regulation (EC) No 1272/2008 [CLP]

Pictogram

Signal word

Danger

#### **Hazard statement(s)**

H318 Causes serious eye damage

#### **Precautionary statements:**

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P313 Get medical advice/attention.

See Section 16 for the full text of H-Statements mentioned in this Section.

#### **2.3 Other hazards**

None known.

### **SECTION 3 COMPOSITION/INFORMATION ON INGREDIENTS**

#### **3.2 Mixtures**

#### **Hazardous ingredients according to Regulation (EC) No 1272/2008**

Component: **Triton X-100** (concentration 10–20%) (included in the Candidate List of Substances of Very High Concern (SVHC) according to Regulation (EC) No 1907/2006 (REACH))  
CAS No: 9002-93-1

EC No: -

A registration number is not available for this substance as the substance or its use are exempted from registration according to Article 2 REACH Regulation (EC) No 1907/2006, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline.

Classification: Acute Tox. 4 (H302); Serious Eye Dam. 1 (H318)

Concentration: 5–10%

See Section 16 for the full text of H-Statements mentioned in this Section.

### **SECTION 4 FIRST AID MEASURES**

#### **4.1 Description of first aid measures**

Consult a physician. Show this safety data sheet to the doctor in attendance.

**If inhaled:** If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

**In case of skin contact:** Take off immediately all contaminated clothing. Rinse skin with water/shower.

**In case of eye contact:** rinse out with plenty of water. Immediately consult an ophthalmologist.

**If swallowed:** immediately make victim drink water (2 glasses at most). Consult a physician.

#### **4.2 Most important symptoms and effects, both acute and delayed**

Irritation and corrosion. Risk of serious damage to eyes.

#### **4.3 Indication of any immediate medical attention and special treatment needed**

No data available.

### **SECTION 5 FIRE FIGHTING MEASURES**

#### **5.1 Extinguishing media**

Suitable extinguishing media: Use water spray, foam, dry chemical or carbon dioxide. (Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.)

Unsuitable extinguishing media: For this substance/mixture no limitations of extinguishing agents are given.

#### **5.2 Special hazards arising from the substance or mixture**

Not combustible. Ambient fire may liberate hazardous vapours.

#### **5.3 Advice for firefighters**

Special protective equipment for firefighters: In the event of fire, wear self-contained breathing apparatus.

Further information: Prevent fire extinguishing water from contaminating surface water or the ground water system.

## TRITON X SOLUTION SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

### **SECTION 6 ACCIDENTAL RELEASE MEASURES**

#### **6.1 Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Do not breathe vapours, aerosols. Avoid substance contact. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see Section 8.

#### **6.2 Environmental precautions**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

#### **6.3 Methods and material for containment and cleaning up**

Soak up with inert absorbent material and dispose of as hazardous waste. Keep in suitable, closed containers for disposal.

#### **6.4 Reference to other sections**

For disposal, see Section 13.

### **SECTION 7 HANDLING AND STORAGE**

#### **7.1 Precautions for safe handling**

Avoid inhalation of vapour or mist. For precautions see Section 2.2.

#### **7.2 Conditions for safe storage, including any incompatibilities**

Tightly closed. Recommended storage temperature see product label.

#### **7.3 Specific end use(s)**

No other specific uses are specified apart from those listed in Section 1.2.

### **SECTION 8 EXPOSURE CONTROLS/PERSONAL PROTECTION**

#### **8.1 Control parameters**

Contains no substances with occupational exposure limit values.

#### **8.2 Exposure controls**

##### **8.2.1 Appropriate engineering controls**

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

##### **8.2.2 Personal protective equipment**

(a) Eye/face protection: Tightly fitting safety goggles. Faceshield (8-inch minimum). Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

(b) Skin protection: Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves should satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

(c) Body Protection: Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

(d) Respiratory protection: Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

#### **8.2.3 Environmental exposure controls**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

### **SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES**

#### **9.1 Information on basic physical and chemical properties**

a) Appearance Form: clear, liquid; Colour: light yellow

b) Odour no data available

c) Odour threshold no data available

d) pH 9.7

e) Melting point / freezing point approx. 6°C

f) Initial boiling point and boiling range 200°C

g) Flash point 251°C

h) Evaporation rate no data available

i) Flammability (solid, gas) no data available

j) Upper/lower flammability or explosive limits no data available

k) Vapour pressure <1 hPa at 25°C

l) Vapour density no data available

m) Relative density 1.070 g/cm<sup>3</sup>

n) Solubility (ies) Soluble in water

o) Partition coefficient: n-octanol/water no data available

p) Auto-ignition temperature no data available

q) Decomposition temperature no data available

r) Viscosity no data available

s) Explosive properties no data available

t) Oxidising properties no data available

#### **9.2 Other information** no data available

### **SECTION 10 STABILITY AND REACTIVITY**

#### **10.1 Reactivity**

No data available.

#### **10.2 Chemical stability**

Stable under recommended storage conditions.

#### **10.3 Possibility of hazardous reactions**

No data available.

#### **10.4 Conditions to avoid**

No data available.

## TRITON X SOLUTION SAFETY DATA SHEET

This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

### **10.5 Incompatible materials**

Strong acids. Strong bases. Strong oxidizing agents.

### **10.6 Hazardous decomposition products**

Other decomposition products—no data available. In the event of fire: see Section 5.

## **SECTION 11 TOXICOLOGICAL INFORMATION**

### **11.1 Information of toxicological effects**

**Acute toxicity:** no data available

**Skin corrosion/irritation:** no data available

**Serious eye damage/eye irritation:** no data available

**Respiratory or skin sensitisation:** no data available

**Germ cell mutagenicity:** no data available

**Carcinogenicity:** IARC: no component of this product present at levels greater than 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

**Reproductive toxicity:** no data available

**Specific target organ toxicity - single exposure:** no data available

**Specific target organ toxicity - repeated exposure:** no data available

**Aspiration hazard:** no data available

**Additional information:** RTECS: not available. To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

### **11.2 Further information**

Triton X-100

Acute oral toxicity: LD50 Rat: 1,800 mg/kg (RTECS)

Germ cell mutagenicity: Genotoxicity in vitro Mutagenicity (mammal cell test): Mouse lymphoma test Result: negative

## **SECTION 12 ECOLOGICAL INFORMATION**

### **12.1 Toxicity**

No data available.

### **12.2 Persistence and degradability**

No data available.

### **12.3 Bioaccumulative potential**

No data available.

### **12.4 Mobility in soil**

No data available.

### **12.5 Results of PBT and vPvB assessment**

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted.

### **12.6 Other adverse effects**

Discharge into the environment must be avoided.

Components: Triton X-100

Toxicity to fish

LC50 Lepomis macrochirus: 2,800 - 3,200 µg/l; 96 h

Toxicity to daphnia and other aquatic invertebrates

LC50 Daphnia magna: 11.2 mg/l; 48 h

### **12.7 Additional information**

No data available.

## **SECTION 13 DISPOSAL CONSIDERATIONS**

### **13.1 Waste treatment methods**

**Product:** Offer surplus and non-recyclable solutions to a licensed disposal company.

**Contaminated packaging:** Dispose of as unused product.

## **SECTION 14 TRANSPORT INFORMATION**

**IATA/DOT/IMDG/TDG:** Not regulated.

**14.1 UN number:** -

**14.2 UN proper shipping name:** -

**14.3 Transport hazard class(es):** -

**14.4 Packing group:** -

**14.5 Environmental hazards:** -

**14.6 Special precautions for user:** -

## **SECTION 15 REGULATORY INFORMATION**

### **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

Substances of very high concern (SVHC): This product does contain substances of very high concern above the respective regulatory limit (>0.1% w/w), Regulation (EC) No 1907/2006 (REACH), Article 57).

Contains: Triton X-100.

### **15.2 Chemical Safety Assessment**

No chemical safety assessment has been carried out for this product.

## **SECTION 16 OTHER INFORMATION**

### **Full text of H-Statements referred to in Sections 2 and 3**

H302 Harmful if swallowed

H318 Causes serious eye damage

Acute Tox. Acute Toxicity

Serious Eye Dam. Serious Eye Damage

Amended sections are indicated by a line in the border.

The information supplied in this SDS is correct to the best of our knowledge. We do not accept any liability for loss, injury or damage, which may result from its use.

# Midi Parasep® SF Faecal Parasite Concentrator

## Code    Product

- |        |   |
|--------|---|
| 149900 | Midi Parasep® SF                              |
| 149910 | Midi Parasep® SF Formalin & Triton X          |
| 149920 | Midi Parasep® SF SAF & Triton X               |
| 149931 | Midi Parasep® SF SafEFix™ Ecological Fixative |
| 249300 | Midi Parasep® SF AlcorFix™                    |

*Discard in accordance with your standard and local operating procedures for clinical waste.*

---

Products can be ordered direct from Apacor or from an appointed distributor

Visit our website for all the latest information [www.apacor.com](http://www.apacor.com) or e-mail: [sales@apacor.com](mailto:sales@apacor.com)

---