

# Serazym® Stool Assay Controls



## Serazym® Stool Assay Controls to be used as external Quality Control of the routine lab

The requirements of quality assurance in the medical laboratory are increasing. To expand the possibilities for monitoring the quality of diagnostic in the laboratory, the concept of the **Serazym® Stool Assay Controls** was developed.



2x fixed shipping dates / year (Spring and Autumn);  
Schedule overview and order form on [www.virotechdiagnostics.com/de/service](http://www.virotechdiagnostics.com/de/service)



Each Control Set consists of 2 vials of stabilized stool suspension for the most important products of the stool diagnostic



Qualitative target values can be found on the homepage of VIROTECH Diagnostics GmbH soon after the 4 week shelflife of the controls

Standing order or individual order within the specified order periods



Production of ordered Control Sets



Shipment via VIROTECH to the Lab



Run and evaluation of the Control Sets during the routine of the Lab.



Publication of qualitative results of target values after expiration of shelflife of the Control Sets (approx. 4 weeks) at [www.virotechdiagnostics.com](http://www.virotechdiagnostics.com) category "service".

### Important:

Control Sets will only be delivered after advance order and are only validated for the use with Serazym® Stool -EIA.

Manufacturer of the Serazym® Product Panel:  
**Seramun Diagnostica GmbH**  
D-15754 Heidesee

**Seramun**  
Diagnostica GmbH

F-Ser-Stool Test Controls-En-02

**VIROTECH Diagnostics GmbH**

📍 Löwenplatz 5 65428 Rüsselsheim am Main

☎ +49 (0)6142 69090

📞 +49 (0)6142 690919

✉ [info@virotechdiagnostics.com](mailto:info@virotechdiagnostics.com)

🌐 [www.virotechdiagnostics.com](http://www.virotechdiagnostics.com)

**VIROTECH**  
DIAGNOSTICS