

Treponema pallidum LINE

(recombinant)

Lues Confirmatory Test in the uniform EcoSystem


Specificity IgG: **99,2%**
Specificity IgM: **98,6%**

Diagnostic Sensitivity,
total assessment
(IgG + IgM): **100%**

Easy and secure handling:

- ▼ Sorted strips in a booklet
- ▼ Clear strip numbering
- ▼ Quality control on each strips
- ▼ Identical incubation times
- ▼ Test results after approx. 100 minutes
- ▼ General components
- ▼ Cut-off control included in the kit

- ▼ Neg./pos. control separately available **-new-**



Order No: WE150G16 or WE150G32, 16 or 32 strips
WE150M16 or We150M32, 16 or 32 strips

F_WE150.00_28032011-13e